

CHALLENGE:

How to create a clean, modern, user friendly clinical setting specifically designed to meet the care-giver’s needs and create an aesthetically pleasing patient environment?

SOLUTION:

The Paladin “HALO” headwall system. With state of the art lighting integration for user specific environmental control and equipment management to facilitate the interface of family/staff/technology.



- Easy to install**
- Single point connections**
- Soft feed gas distribution to outlets**
- Used with Fairfield-style equipment rail**



- **Focused on creating the ideal clinical setting**
- **Flexibility to modify clinical needs while preserving the value of the wall system**
- **Adaptability to cost effectively manage changing technology and user needs**
- **Expandability to quickly add essential services during surge and other unforeseen events**
- **Independent and flexible state of the art lighting integration**
- Custom configuration of gas and electrical outlets
- Specify any manufacturer or type of outlet
- Any Marlite FRP finish for front panel
- Non-invasive gas expansion with no system certification required
- Integrated with “Evolution” Fairfield-style equipment rail system for ideal equipment management to clear the floor of technology and rolling stands. Improves patient access, reduces trip hazards, and improves floor cleaning access.
- LED color changing technology for choice in lighting color spectrums