



CHALLENGE:

How to create a clean, modern, user friendly environment with a cost effective and simple to install headwall product?

**SOLUTION:**

The Paladin “Revolution” headwall system is a state of the art integration of lighting, infection control, equipment management, and simple installation technology.

**Easy to install****Single point connections****Soft feed gas distribution to outlets****Double layers of integrated “EVOLUTION” equipment rail**

- **Flexibility to modify clinical needs while preserving the value of the wall system**
- **Adaptability to cost effectively manage changing technology and user needs**
- **Expandability to quickly add essential services during surge and other unforeseen events**
- **Integrated with “Evolution” Fairfield-style equipment rail system. Equipment rail clears the floor of technology and rolling stands which improves patient access, reduces trip hazards, and improves floor cleaning access.**
- Custom configuration of gas and electrical outlets
- Specify any manufacturer or type of outlet
- White with embedded anti-microbial finish, or a wide variety of RAL paint or wood grain finishes
- Supports hundreds of equipment management tools and accessories
- Non-invasive gas expansion with no system certification required