



	ITEM CELL	ITEM DESCRIPTION Item numbers are dependent upon exact product specifications and may not be listed below.
 <p>VARIABLE AND NOT SHOWN ON RAIL FOR APPROXIMATE PLACEMENT ON RAIL..</p> <p>GCX CHANNELS CAN BE MOUNTED TO 2 OR MORE RAILS.</p> <p>THERE MUST BE AN OVERLAP OF APPROXIMATELY 2-4" ABOVE TO RAIL AND BELOW BOTTOM RAIL.</p>		<p>GCX149 Rail Mount Kit for GCX style Swing Arms Integrates small 4" GCX channel to the Equipment Rail Pre-fabricated to include GCX channel and rail mounting adapter Does not include patient monitor or swing arm mount The ONLY single rail mount option for GCX available. All other solutions use 2 or more rails.</p> <p>GCX CHANNEL For mounting various equipment such as: computer monitors, keyboards, CPU units, patient monitors, lighting, IV's, etc.</p> <p>Length of channel and number of rails used determines the Paladin product number.</p>