






















<p>SHOWN ON RAIL FOR APPROXIMATE PLACEMENT POINT ON RAIL.</p> <p>VACUUM BOTTLE PLACEMENT IS VARIABLE DEPENDING ON PRODUCT ITEM AND VENDOR (3) EXAMPLES SHOWN</p>	<p>ITEM CELL</p>	<p>ITEM DESCRIPTION</p>
         	          	<p>AC101-AM Vacuum bottle slide, Universal for all vacuum bottle sizes/manufacturers anti-microbial finish</p> <p>AC101-AM Vacuum bottle slide, Universal for all vacuum bottle sizes/manufacturers anti-microbial finish</p> <p>AC101-AM Vacuum bottle slide, Universal for all vacuum bottle sizes/manufacturers anti-microbial finish</p> <p>GBA-101 Duplex Gas Block for Medical Air for splitting Air outlet shown with flowmeters, but does not include flowmeters</p> <p>GBO-101 Duplex Gas Block for Medical Oxygen for splitting Oxygen outlet shown with flowmeters, but does not include flowmeters</p> <p>Vacuum Suction Regulator Single vac suction regulator does not often directly rail mount Normally it mounts directly in to the outlet unless the outlet needs to be split to double capacity. Then see below.</p> <p>GBV-101 Vacuum Suction Regulator Single vac suction regulator does not often directly rail mount Normally it mounts directly in to the outlet unless the outlet needs to be split to double capacity. Then see below.</p> <p>AC104-AM Blender bracket, universal for any brand of blender Does not include Blender Does not show tubing or other associated equipment</p> <p>VM201 Double ventilator manifold with DISS outlets, 2 Air and 2 Oxygen Does not show or include the Blender that can be mounted into the attached AC101 on the face of the unit</p> <p>AC105 Cylinder cage and mount for D or E tank cylinders shown with Medical Air, does not include tank</p> <p>AC105 Cylinder cage and mount for D or E tank cylinders shown with Medical Oxygen, does not include tank</p>