

PACKAGE:
----------

APPROVAL DATE:

# PALADIN HEALTHCARE RAIL SYSTEM TYPICAL SOLUTION

SCALE: N.T.S.

REV 2:

DATE: 3-24-20


REV 3:

REV 1: 4-13-20

REV 4:

## IE.2

SHEET NO.



Item	Description	Qty
R204	4 Ft. (48") Equipment Rail, pre-drilled, Including End Caps and Laminate Inlay	1
GBA105	Air Gas Block (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
GBO105	Oxygen Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
GBV105	Vacuum Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
GHO201-3-0	Oxygen Hose, 3 ft Length, DISS Female Hand Nut x DISS Female Hand Nut	1
GHA201-3-0	Air Hose, 3 ft Length, Diss Female to Diss Female, straight, (36", Clear, Non-Conductive)	1
GHV201-3-0	Vacuum Hose, 3 ft Length, Diss Female to Diss Female Handnut, Straight (36", Clear, Non-Conductive)	1
AC101-ABS-Fair	Vacuum Bottle Slide, incl. rail mount kit, White, Anti-Microbial ABS with Anti-Microbial Adapter	2



940 W Oakland Avenue, Suite A-1  
Oakland, FL 34787  
407.614.8781 Office / 407.614.8782 Fax

**ELEVATION: 4ft Gas Solution**

**REF CODE: QB4522**

Description	Qty
4 Ft. (48") Equipment Rail, pre-drilled, Including End Caps and Laminate Inlay	1
Air Gas Block (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
Oxygen Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
Vacuum Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	1
Oxygen Hose, 3 ft Length, DISS Female Hand Nut x DISS Female Hand Nut	1
Air Hose, 3 ft Length, Diss Female to Diss Female, straight, (36", Clear, Non-Conductive)	1
Vacuum Hose, 3 ft Length, Diss Female to Diss Female Handnut, Straight (36", Clear, Non-Conductive)	1
Vacuum Bottle Slide, incl. rail mount kit, White, Anti-Microbial ABS with Anti-Microbial Adapter	2
For a formal quotation specific to your project, please contact us at: admin@paladinhc.com or 407.614.8781 Office or 407.614.8782 Fax	