



Item	Description	Qty
R206	6 Ft. (72") Equipment Rail, pre-drilled, Including End Caps and Laminate Inlay	1
GBA105	Air Gas Block (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
GBO105	Oxygen Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
GBV105	Vacuum Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
GHO201-3-0	Oxygen Hose, 3 ft Length, DISS Female Hand Nut x DISS Female Hand Nut	2
GHA201-3-0	Air Hose, 3 ft Length, Diss Female to Diss Female, straight, (36", Clear, Non-Conductive)	2
GHV201-3-0	Vacuum Hose, 3 ft Length, Diss Female to Diss Female Handnut, Straight (36", Clear, Non-Conductive)	2
SBA201	Air Outlet Splitter Block DISS Female Quick connect to (2) DISS Male Outlets with C/V	1
SBO201	Oxygen Outlet Splitter Block DISS connection to (2) DISS Outlets	1
SBV201	Vacuum Outlet Splitter Block DISS Quick connect to (2) DISS Outlets	1
AC101-ABS-Fair	Vacuum Bottle Slide, incl. rail mount kit, White, Anti-Microbial ABS with Anti-Microbial Adapter	4



**PALADIN
HEALTHCARE**

940 W. Oakland Avenue, A1
Oakland, FL 34787
(407) 614-8781 Office
gschindele@paladinhc.com

PACKAGE:

APPROVAL DATE:

PALADIN HEALTHCARE RAIL SYSTEM TYPICAL SOLUTION

SCALE: N.T.S.

REV 2:

DATE: 3-24-20

REV 3:

REV 1: 4-13-20

REV 4:

IE.1

SHEET NO.

C:\Users\lra\OneDrive\Paladin Healthcare\Projects 2020\NAVY\Paladin Typical Solution NAVY_Rev1.dwg, 4/14/2020 9:59:20 AM, DWG To PDF.pc3

COPYRIGHT © 2013 BY PALADIN HEALTHCARE LLC. ALL RIGHTS RESERVED. THIS DOCUMENT IS THE PROPERTY OF PALADIN HEALTHCARE LLC. IT IS TO BE USED ONLY FOR THE SPECIFIC PROJECT IDENTIFIED OR REFERRED TO HEREIN AND IS NOT TO BE USED ON OTHER PROJECTS OR EXTENSIONS TO THIS PROJECT, IN WHOLE OR PART, EXCEPT BY EXPRESS WRITTEN AGREEMENT WITH AND COMPENSATION TO PALADIN HEALTHCARE LLC.



940 W Oakland Avenue, Suite A-1
Oakland, FL 34787
407.614.8781 Office / 407.614.8782 Fax

ELEVATION: Typical 6ft Gas Solution

REF CODE: QB 4520

Description	Qty
6 Ft. (72") Equipment Rail, pre-drilled, Including End Caps and Laminate Inlay	1
Air Gas Block (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
Oxygen Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
Vacuum Gas Block, (2) 1/8" NPT front outlet, 90 degree DISS (CV) center outlet	2
Oxygen Hose, 3 ft Length, DISS Female Hand Nut x DISS Female Hand Nut	2
Air Hose, 3 ft Length, Diss Female to Diss Female, straight, (36", Clear, Non-Conductive)	2
Vacuum Hose, 3 ft Length, Diss Female to Diss Female Handnut, Straight (36", Clear, Non-Conductive)	2
Air Outlet Splitter Block Diss Female Quick connect to (2) DISS Male Outlets with C/V	1
Oxygen Outlet Splitter Block DISS connection to (2) DISS Outlets	1
Vacuum Outlet Splitter Block DISS Quick connect to (2) DISS Outlets	1
Vacuum Bottle Slide, incl. rail mount kit, White, Anti-Microbial ABS with Anti-Microbial Adapter	4
For a formal quotation specific to your project, please contact us at: admin@paladinhc.com or 407.614.8781 Office or 407.614.8782 Fax	