



**Rail Mounting Kit with  
Medical Chart Holder,  
HIPAA compatible, vertical,  
Available in clear front as well**

**Item Numbers:  
AC174-2-PM-HIPAA (2-tiers)  
JSN Number: not established**

**Dimensions:  
12.5"W x 12.375"H x 3.5"D  
Weight: 5lbs**

**Item Description:**

**Photo shows 3-tiers chart holder for reference**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.