



** Welch Allyn
equipment SOLD SEPARATELY

Combination Mounting kit

for both:

Owner provided
Welch Allyn Ophthalmoscope

(models 767/777) and

Owner provided
Specula holders

Item Number: DEM198

JSN Number: M4200-RM

Dimensions Specula Rail Mount Kit:
2.5”H x 6”W x 0.125”D - Weight: 1lb

Dimensions Ophthalmoscope
Rail Mount Kit:

6”H x 2.5”W x 0.125”D - Weight: 1lb

Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.