



**Universal Face Mask
Dispenser Holder, including rail
mounting kit, heavy duty
acrylic**

**For a complete overview which
brand items will fit, please see
next page or visit
www.PaladinHC.com**



**Item Number: IC-PFP021-PM
JSN Number: not established**

**Dimensions:
8.125"W x 7.5"H x 5.625"D
Weight: 3lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.



PART NUMBER: IC-PFP021-PM – IC219PM

DESCRIPTION: Universal Face Mask Dispenser Holder

NOTES: 8.125"W x 5.625"D x 7.5"H (dispenser only)
Heavy Duty Acrylic, Dispenser only,
incl. Rail Mount Kit.

**FITS THE
FOLLOWING
DISPENSERS:**

Cardinal Health: AT7511, AT71021, AT771141, AT771211, AT771212,
N95-ML

Halyard: 25867, 25868, 25869, 28797, 28806, 28808, 28810, 32856, 37525,
41802, 41805, 47080, 47085, 47090, 47107, 47108, 47117, 47119, 47127,
47295, 47297, 48100, 48105, 48201, 48215, 48296, 49214, 49310, 49700,
49701, 59927, 59928, 59929, 59931, 59936, 59938, 59939, 59940, 62113,
62115, 62356, 62357, 62362, 62363, 62364

McKesson: 91-1200, 91-2002

MediChoice, Primagard/PriMed: PG41001, PG41092, PG41200, PG41201,
PG41241, PG41331, PG42092, PG42121, PG42332, 205ELM001

Medline: NONSH400, NON27110, NON27120, NON27120-2, NON27122,
NON27355, NON27365, NON27374, NON27376, NON27377, NON27378,
NON27385, NON27402, NON27408EL, NON27412, NON27412EL,
NON27600, NON27600AG, NON27710, NON27710EL, NON27712EL,
NON27810, NON27810EL, NON27812EL,

Precept: 65-3320, 65-3322, 65-3344, 15100, 15101, 15130, 15150, 15155,
15200, 15211, 15212, 15230, 15240, 15260, 15300, 15310, 15320, 15330,
15420, 15500, 15710

Prestige Ameritech/ProGear: FM81000, FM82120, FM84210, FM85011

Sperian: 14110452

Tidi: 9010, 9020, 9030, 9040