



**Medical Gas Block with**  
**2 x Female Quick Connect front**  
**outlets. Available with Ohio,**  
**Chemetron, Oxyquip and Schraeder**  
**fittings.**

**Air and Oxygen blocks come with an additional 90 degrees DISS center outlet check valve. Vacuum block does not include center check valve.**

**Item Numbers:**

**GBA103 (Air)**

**GBO103 (Oxygen)**

**GBV103 (Vacuum)**

**JSN Number: not established**

**Dimensions: 4"W x 3"H x 3.5"D**

**Weight: 3lbs**

**Flow meters and vacuum regulators are not included**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.