



PART NUMBER: HWH-006 / HWH-008 / HWH-009 (JSN # not established)

DESCRIPTION: REVOLUTION Horizontal Headwall Unit

- NOTES:**
- Available in 6', 8', and 9' lengths, or customized to project specs
 - Customizable Gas and Electrical outlet type and placement
 - (2) layers of integrated "Evolution" Fairfield-style equipment rail
 - Single-point connection feed manifold
 - Option for integrated lighting, up and/or down
 - Optional RAL paint finishes available for face panels
 - Any color Wilsonart laminate inlay standard for rails. Alternate manufacturers available.
 - Fully integrates with Marlite wall protection system

