



**PART NUMBER:** HWH-006 / HWH-008 / HWH-009 (JSN # not established)

**DESCRIPTION:** REVOLUTION Horizontal Headwall Unit

- NOTES:**
- Available in 6', 8', and 9' lengths, or customized to project specs
  - Customizable Gas and Electrical outlet type and placement
  - (2) layers of integrated "Evolution" Fairfield-style equipment rail
  - Single-point connection feed manifold
  - Option for integrated lighting, up and/or down
  - Optional RAL paint finishes available for face panels
  - Any color Wilsonart laminate inlay standard for rails. Alternate manufacturers available.
  - Fully integrates with Marlite wall protection system

