



**Rail Mounting Kit for owner
provided Daniels “Sharpsmart”
Container S14 + S22**

**Item Number:
IC304-P (locking)
IC309-P (non-locking)
JSN Number: A5108**

**Dimensions:
4”W x 0.125”D x 3.375”H (locking)
7”W x 0.125”D x 4”H (non-locking)
Weight: 2lbs**

**** Does not include sharps container
Specific sharps container (and locking box)
model number(s) must be verified prior to order
placement**

Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.