



Rail Mounting Kit with 4 bins,
top feed, surgical PPE
dispenser, heavy duty acrylic,
fits bulk face masks, shoe
covers, etc. bottom of
dispenser slotted for easy
cleaning

Item Number: IC543-PM

JSN Number: not established

Dimensions:

23.75"W x 7"H x 5.5"D

Weight: 5lbs

*Customer to define mounting preference
(horizontal or vertical) at time of order.*

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.