



Medical-Gas Block
with 2 x 1/8" NPT Male front
outlets. Air and Oxygen blocks
come with an additional 90
degrees DISS center outlet check
valve. Vacuum block does not
include center check valve.

Item Numbers:

GBA101 (Air)

GBO101 (Oxygen)

GBV101 (Vacuum)

JSN Number: not established

**Flow meters and vacuum regulators
are not included**

Dimensions: 4"W x 3"H x 3.5"D

Weight: 3lbs

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.