



Rail Mounting Kit

for owner provided FERNO bags.

Available in 3 different colors:

Red, white and blue

Item Numbers:

SM501-P (Plate only)

PR-FN504B (medium bag)

PR-FN504R (large bag)

JSN Number: not established

Dimensions:

11”H x 13”W x 1”D (plate only)

13”H x 18”W x 9”D (small bag)

19”H x 18”W x 9”D (large bag)

Weight: 2lbs (plate only)

5lbs (small bag), 6lbs (large bag)



Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.