



MADE IN THE USA


Architectural Solutions for Clinical Environments with the Fairfield-Style Rail



PALADIN HEALTHCARE®

CHALLENGE

How to Future-Proof Healthcare Spaces for the best Patient Care and Caregiver Safety?



SOLUTION

Equipment Rail can be your insurance against future changes. Paladin Healthcare LLC is your value-add partner, bringing over 50 years of experience in this field to your project. **If it's on a wall – It's on a rail!**

EVOLUTION -

Equipment Management Rail (Fairfield-style)

- Flexible, Adaptable, Expandable Equipment and Materials Solutions for Life-cycle Cost Reductions
- Easiest to use adapter mounts which may be locked to the rail for security
- Embedded Anti-microbial features in adapters and plates
- The largest readily available selection of laminate inlay colors in the industry
- Made in the USA

REVOLUTION

Horizontal Headwall System

- 2 layers of "Evolution" Fairfield-style rail extruded onto face frame
- Any manufacturers outlets for gas and electrical supply
- Standard or custom placement of all services
- Single manifold with 3-point brazing connections
- Anti-microbial powder coat finish: white or grey (standard), or any RAL paint color
- Made in the USA

REVOLUTION - V

Vertical Headwall System

- Designed to support easy clinical access and high-style aesthetics
- Ergonomic gas and electrical services on the face of the headwall
- Any manufacturers outlets for gas and electrical supply
- Standard or custom placement of all services
- Single manifold with 3-point brazing connections
- Removable faceplate for easy service and future decor changes
- Made in the USA

F.T.S.

Forward Treatment Station and Fast Transport Cart

- Forward Triage Cart is the ultimate Emergency multi-purpose tool for Forward Triage in Waiting Rooms, Mass Casualty/Surge Patient Management, and Psych Rooms
- Fast Transport Cart attaches to the stretcher for management of all equipment and supplies during patient transport
- Both models are fully customizable with hundreds of management tools and accessories, including portable monitors and power cables
- Made in the USA





About us...



Paladin Healthcare LLC is now a Vizient Supplier!
Vizient is the industry's largest supply contracting company and the nation's largest member-owned health care services company.
Are you a Vizient member yet? Vizient Supplier Number: **EC7040**



Paladin Healthcare LLC is now a Premier Supplier!
With the network of 4.000+ Hospitals and Health Systems you will now be able to add this Premier Agreement to your Sales Toolkit!
Paladin-Premier-Contract Number: **CO-FA108**



You can find Paladin Healthcare LLC on www.attainia.com.



Paladin Healthcare LLC is an Ascension Supplier!
Ascension provides care in over 19 states and the District of Columbia.
Ascension Contract Number: **0000090255**



Paladin Healthcare LLC is now an Advent Health Supplier! Advent Health provides care in 9 states. Advent Health Contract Number: **AH-FA-264**



You can find Paladin Healthcare LLC on www.dla.mil. ECAT is an online ordering, distribution and payment system providing DoD and other federal customers access to multiple manufacturer and distributor commercial catalogs at discounted prices. DLA Troop Support's Medical supply chain developed ECAT to streamline its business practices and expand its range of procurement options.



Paladin Healthcare LLC offers an education for your team to learn more about "What happens after the open house? - What can we do to manage clinical reality?" and "Mass Casualty – How to plan for the unknown"
Our Lunch & Learns can take place in your offices or in our brand-new state of the art learning center in Oakland, FL (Orlando).
These courses are registered for 1 hour of continuing education credit each with the American Institute of Architects (AIA) and The Center for Health Design (EDAC).



You can find a complete library of CAD files and a growing number of BIM / REVIT files through our members log-in section on www.paladinhc.com.