



PART NUMBER: HWV-HALO (different models available) (JSN # not established)

DESCRIPTION: HALO Vertical Headwall Unit

NOTES:

- Vertical headwall unit that can be used alone as a headboard for a single bed or for use with two beds
- Ergonomic gas and electrical services to promote clinical excellence
- Single-point connection manifold. Only 3 brazing connections.
- Customizable width and services placement available
- Standard bundled solution includes (6) 72" "Evolution" equipment rails
- Any color Wilsonart laminate inlay standard for rails. Alternate manufacturers available.
- Face panel customizable in Marlite FRP finishes
- Fully integrates with Marlite wall protection system

