



# “HALO” Headwall System



## CHALLENGE:

How to create a clean, modern, user friendly clinical setting specifically designed to meet the caregiver’s needs and create an aesthetically pleasing patient environment?

## SOLUTION:

The Paladin “HALO” headwall system. With state of the art lighting integration for user specific environmental control and equipment management to facilitate the interface of family/staff/technology.



- Easy to install**
- Single point connections**
- Soft feed gas distribution to outlets**
- Used with Fairfield-style equipment rail**



- **Focused on creating the ideal clinical setting**
- **Flexibility to modify clinical needs while preserving the value of the wall system**
- **Adaptability to cost effectively manage changing technology and user needs**
- **Expandability to quickly add essential services during surge and other unforeseen events**
- **Independent and flexible state of the art lighting integration**
- Custom configuration of gas and electrical outlets
- Specify any manufacturer or type of outlet
- Any Marlite FRP finish for front panel
- Non-invasive gas expansion with no system certification required
- Integrated with “Evolution” Fairfield-style equipment rail system for ideal equipment management to clear the floor of technology and rolling stands. Improves patient access, reduces trip hazards, and improves floor cleaning access.
- LED color changing technology for choice in lighting color spectrums





**“HALO”**  
with **“EVOLUTION”** Equipment Rail  
for Lifecycle Cost Management



**PALADIN**  
**HEALTHCARE®**

**“HALO” WITH  
“EVOLUTION”**

Promote lifecycle flexibility with the combined headwall and rail management systems.

Equipment, accessories, and materials will be changed, upgraded, and relocated at will and as required.

Gases and Electrical placement are standard or customized.

Face panel can be easily changed to meet aesthetic changes in future renovations.

State of the art lighting is fully user controlled.

Further promote clinical excellence and infection control with the use of “Evolution” on the side walls for maximized materials management. Place things:

- Where They Need It
- When They Need It
- So They Can Easily Find It

Permanently end unnecessary wall damage due to changing infection control product contracts.



“HALO” vertical headwall with removable faceplate to expose gas and electrical services for quick and easy access. Customizable outlet types and locations. Integrated LED lighting and optional nightlight at base. Can be made as a chase to discretely plumb services from the ceiling for an attractive solution to improved construction costs, simplified coordination, and lifetime maintenance.

