


Architectural Solutions for Clinical Environments



PRODUCT CATALOG

10-20-2023



Accessories available for Modular, Hill-Rom/Wittrock and AMICO Style Rails



Table of Contents (Online)

Accessories for Modular, Hill-Rom/Wittrock, AMICO	2
Evolution – Fairfield Style Equipment Management Rail	6
Equipment Management	10
Infection Control Devices	17
EVS Environmental Services	19
Soap / Sanitizer / Wipes	21
Paper Towel Dispensers	25
Glove Box Dispensers	28
O.R. / Scrub Sink Solutions / PPE Dispensers	31
Sharps Containers	42
Baskets and Bins	44
Clips, Hooks and Wire Management	49
Shelves and Work Surfaces	53
Monitor Mounting and Cable Management Tools	56
Medical Gas and Vacuum Management	62
Infusion Management and Lighting	69
Miscellaneous Accessories	72
Architectural Headwall Systems	80
Mobile Treatment and Transport Systems	86
Critical Care Anywhere	92
Equipment Management	93
Optimize your Patient Care Environment	94
Transform 40 year old ED Bays	95
Lean Process Improvement	96
Handwashing, HAI's and Reimbursement	97
True Terminal Cleaning and Infection Control	98
Complete Design Services	99

Architectual Solutions for Clinical Environments



Hospital Grade
Power Strip



Holder for round
bottle type dispensers



Rail Mount for First Aid Kits



Holder for
Crutches



Rail Mounts for Artwork



Rail Mount for
Owner provided
Scrub Holder



Rail Mount for
different brands
hand sanitizers



Rail Mount
for
Thermometer



Rail Mount
for Mirror



How to read our part numbers (*very important*)

Part numbers ending in “-P” = Mounting Plate (Rail Mounting Solution without Accessory)

This means the item to be mounted is NOT included in the purchase. It's just the mounting plate and adapter kit necessary to mount the item to the rail. Some additional installation maybe required to a 3rd party device.

Part numbers ending in “-A” = Adapter (Rail Mounting Solution without Accessory)

This means the item to be mounted is NOT included in the purchase. Just the adapter (no plate) will be delivered. The adapter can be mounted directly on the accessory. Some additional installation maybe required to a 3rd party device.

Part numbers ending in “-PM” = Product Mounted (to Mounting Plate)

This means the product itself is Paladin factory mounted to the plate and adapter kit necessary to mount the item to the rail. The item and adapter kit are both included in the price. Maybe used for its intended purpose out of the box.

Part numbers including “HV” = Horizontal or Vertical (mounting options)

This means the mounting plate has integral options for mounting the item either horizontally, or vertically, on the rail.

Part numbers ending in “-K” = Kit (2 or more Rail Mounting Solution in a Kit)

This means it is a top-level part number incorporating several mounting plates with or without accessories attached (see item description for more details) to fulfill specific requirements.

How our Evolution rail system works

For an image-based explanation, see page 6: “Required Parts and Pieces”

We use a 3-Part system to turn your wall mounted items into a componentized method of mounting and managing your equipment in a way that is flexible, adaptable, and expandable now and throughout the lifecycle of the facility.

1. Rail (Evolution Equipment Rail)
2. Plate + Adapter = Mounting Kit
3. traditionally wall mounted Item

The mounting kit marries the traditionally wall-mounted item to the rail. Items can be locked in place on the rail by tightening the Allen screws on the bottom of the adapter.

When a change of devices or process is necessary, simply snap the item off of the rail, snap it onto the rail in its new location, test your new workflow process, and lock the items to the rail when the layout is efficient. Change is a Snap!

We have hundreds of product mounting solutions available

If you don't find what you need in this catalog, please call us! We might have a solution already! If we haven't, we are more than happy to design one for you!

Evolution



Fairfield-Style

Medical Equipment Management Rail and Accessories System




EVOLUTION: Required Parts and Pieces




Paladin Healthcare's Fairfield-style Equipment Rail System



Adapters – every “mounting kit” includes an adapter





****The Mounting Kit is determined by the device to be hung on the rail. User provides device specification.****





Most Flexible Product Available in the Industry



Paladin Healthcare unique and “salient” features are as follows:

1. Anti-Microbial finish on the plates, adapters and wall stand-offs.
2. Wall stand-off to reduce leaning and the “bounce” effect when using a device.
3. Multiple holes to attach future products from different vendors without the need to purchase new plates each time.
4. The ability to adjust the rail mounting adapter to different heights relative to the fixed rail, allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.



Pictures below show:

1 item – same adapter - 4 different heights





MADE IN THE USA




FOR OFFICIAL CSI MASTER FORMAT SPECIFICATION AND DETAILED INSTALLATION INSTRUCTIONS, PLEASE VISIT OUR WEBSITE.

SPECIFICATION: EVOLUTION Fairfield-Style Equipment Management Rail PART NUMBER: R-2XX


SPECIFICATION: SOLUTION Universal Equipment Management Rail PART NUMBER: MR-2XX



** Image shown is a seamless fit in a modular casework application and does not include end caps. **



Evolution Fairfield Style Rail



Solution Universal Rail

- Paladin Equipment Management Rails: a seamless, extruded aluminum profile with clear anodized finish containing a groove for attaching an adapter that supports the Owner supplied specified hospital accessories and equipment and includes termination caps at each end. The groove and adapters are designed so that the adaptor requires a seven-degree tilt for attachment. The rail is designed to accommodate a seamless Wilsonart laminate decorative color inlay the length of the rail up to ten feet in length.
- Paladin Equipment Management Rails shall be specified in one to ten-foot lengths as shown on drawings or indicated on purchase order to accommodate room conditions. Rail will be cut to exact lengths as indicated. For lengths greater than ten feet, a JOINER piece is used to attach multiple pieces of rail to meet the required total length. **The length of the rail is nominal for the rail only and does *not* include the end caps. Use of two end caps increases the total length of the rail by ½”.
- Paladin Equipment Management Rail is pre-drilled on sixteen-inch centers (minimum) with the exception of a one-foot rail length which has two mounting holes drilled at eight inches on center. The Rail may also be field drilled by the installer to accommodate specific job site conditions.
- End caps shall be constructed of injected molded thermoplastic and shall resist most hospital cleaning and disinfecting chemicals. The color shall be gray and molded throughout the material. End caps shall be friction-fit without use of fastening devices.
- Decorative inlay shall be supplied by Paladin for installation by the owner or contractor. Inlay is Wilsonart laminate (www.Wilsonart.com) in the standard finish which will resist most hospital cleaning and disinfecting chemicals. Color shall be specified by the Owner. Specialty coating finishes can be made available for additional cost and lead time. Alternate laminate manufacturers may be specified for additional cost and lead time. Contact manufacturer for details.

Equipment Management

Rail Mounting Solutions

For Diagnostic Instruments - Any Make and Model

Some Devices may also be available upon request





MADE IN THE USA


PART NUMBER: DEM106-P
(1.5"H x 5.5"W x 1"D)

DESCRIPTION: Mounting Plate for owner provided Welch Allyn / Hillrom aneroid BP sphygmomanometer Tycos 509 (Welch Allyn / Hillrom product is not included, shown with cuff basket SM101-PM not included)



PART NUMBER: DEM116-A
(3.5"H x 3.5"W x 1"D)

DESCRIPTION: Adapter for owner provided Baum Aneroid BP sphygmomanometer (Baum product is not included)



PART NUMBER: SM101-PM
(7"H x 6.75"W x 4"D)

DESCRIPTION: Mounting Plate with BP cuff basket factory attached

NOTES: Plastic coated wire U-shaped basket. Max load 4lbs.



PART NUMBER: DEM103-PM
(1.5"H x 7.5"W x 4.75"D)

DESCRIPTION: Mounting Plate with Welch Allyn / Hillrom 71420 tray for owner provided 71140 desk charger (charger and instruments are not included)






MADE IN THE USA

PART NUMBER: DEM198-K

DESCRIPTION: Combination Mounting Kit:
Mounting plates for both owner
provide Welch Allyn / Hillrom
ophthalmoscope (models
767/777) and owner provided
specula holders (Welch Allyn /
Hillrom products are not
included)




PART NUMBER: DEM108-P
(6"W x 2.5"H x 1"D)

DESCRIPTION: Mounting Plate for owner
provided Welch Allyn / Hillrom
ophthalmoscope (models
767/777) (Welch Allyn / Hillrom
products are not included)



PART NUMBER: DEM109-P
(2.75"W x 4"H x 1"D)

DESCRIPTION: Mounting Plate for owner
provided Welch Allyn / Hillrom
specula holders (Welch Allyn /
Hillrom product is not included)



PART NUMBER: DEM110-P
(7.5"W x 4.25"H x 1"D)

DESCRIPTION: Mounting Plate for owner
provided Welch Allyn / Hillrom
"Sure Temp" thermometer
model 692 (Welch Allyn /
Hillrom product is not
included)





MADE IN THE USA

PART NUMBER: DEM170-P
(1.5"W x 6"H x 1"D)

DESCRIPTION: Mounting Plate for owner provided wall mount for Welch Allyn / Hillrom thermometer Thermoscan Pro 6000, or for Genius 3 Tympanic thermometer #CESS-A-GEN030 (Welch Allyn / Hillrom products are not included)



PART NUMBER: DEM146-A
(3"W x 2.6"H x 1"D)

DESCRIPTION: Adapter for owner provided Exergen Temporal Scanner Thermometer TAT5000 or wall mount Exergen 134201



PART NUMBER: DEM133-P
(2"W x 5.5"H x 1"D)

DESCRIPTION: Mounting Plate for owner provided Exergen Wall Mount 134305 and 134306



PART NUMBER: DEM100-K-4

DESCRIPTION: Complete Mounting Plate Kit with 4ft rail, adapters and BP cuff basket for the equipment components of the owner provided Welch Allyn / Hillrom 777 diagnostic instruments on a 4ft Rail (Welch Allyn / Hillrom products are not included)





MADE IN THE USA

PART NUMBER: DEM196-K

DESCRIPTION: Combination Mounting Kit includes Mounting Plates for both owner provided Welch Allyn / Hillrom ProBP 3400 Model wall mount #34XFWT-B and BP cuff basket factory attached (Welch Allyn / Hillrom product is not included)



PART NUMBER: DEM105-P
(7"H x 11"W x 1"D)

DESCRIPTION: Mounting Plate incl. wall stand offs for owner provided Welch Allyn / Hillrom "Connex Integrated Wall System" comprehensive diagnostic center (Welch Allyn / Hillrom products are not included)



PART NUMBER: DEM164-P
(3"W x 12"H x 1"D)

DESCRIPTION: Mounting Plate for owner provided Welch Allyn / Hillrom S-shaped spot monitors with horizontal monitor mount (Welch Allyn / Hillrom products are not included)



NOTES: Various part numbers available for this monitor type
Part Number: 42NTB-E1 - GTIN: 00732094005714

PART NUMBER: DEM220-K

DESCRIPTION: Mounting Plate Kit for owner provided Kendall SCD700 series

NOTES: The Kit includes SM309 heavy duty U-shape hook, SM306 Cord Cleat and SM301-2 Single Hook 2"





MADE IN THE USA

PART NUMBER: DEM167-PM
(4"W x 6"H x 8"D)

DESCRIPTION: Mounting Plate with monitor arm and GCX Channel
factory mounted for owner provided Welch Allyn / Hillrom 300 series vital signs monitor

(Welch Allyn / Hillrom products are not included)



PART NUMBER: DEM192-K

DESCRIPTION: Combination Mounting Kit includes Mounting Plates for owner provided Welch Allyn / Hillrom Spot Monitor Model 74C and Model 4400 charging pack. The kit includes:
1x mounting plate for spot monitor
1x mounting plate for charging pack



PART NUMBER: DEM161-P
(4"W x 4.5"H x 1"D)

DESCRIPTION: Universal Mounting Plate for owner provided Bair Paws warming unit and the owner provided 3M Bair Hugger warming unit.





MADE IN THE USA

PART NUMBER: DEM210-P
(3.5"W x 6"H x 1"D)

DESCRIPTION: Mounting Plate
for owner provided Midmark
HEINE B-000-11-149-166
AllSpec Disposable Otoscope
Tip Dispenser



PART NUMBER: DEM211-P
(8"W x 4"H x 1"D)

DESCRIPTION: Mounting Plate
for owner provided
Midmark EN200
HEINE Diagnostic Wall Kit



PART NUMBER: DEM212-A
(3"W x 2.6"H x 1"D)

DESCRIPTION: Adapter
for owner provided Midmark
HEINE GAMMA XXL LF
Sphygmomanometer



PART NUMBER: DEM181-P
(6"W x 6"H x 1"D)

DESCRIPTION: Mounting Plate for owner
provided Riester vital signs monitor



Infection Control Devices

Rail Mounting Solutions

EVS Environmental Services

Soap / Sanitizer / Wipes

Paper Towels

Gloves / PPE




Sharps Containers





MADE IN THE USA

*Published by The
Center for Health
Design 2011*



Harley Ellis Devereaux

Standard Process Handwashing System in Patient Rooms;
All new patient care units

EBD Goal:

To create a standard process that facilitates the practice of appropriate hand washing prior to caregiver interaction with the patient.

PROBLEM: Per the CDC, many studies have shown that **the bacteria that cause hospital-acquired infections are most frequently spread from one patient to another on the hands of healthcare workers.**


In 34 studies of hand washing, workers washed their hands only 40% of the time.

SOLUTION: Create a standard directional flow with easy access to sink and hand washing accessories at the room entry. **“...placing the hand washing accessories on an equipment rail allows adjustment in the system without creating damage to the wall and potential home for pathogens.”**

Rail Mounting Solutions

Paladin Healthcare has created a flexible solution for Environmental Services rooms and closets. Our rail and accessories help create a safe and ergonomically appropriate storage solution for your EVS teams. Keep items off the floor, and work surfaces free of clutter by managing your EVS tools (including disinfection solution mixing cabinets) on the rail without redundant damage to the walls.

If it's on a wall, it's on a rail....even in your environment services spaces.





MADE IN THE USA

PART NUMBER: AC285-1-PM Single
(2.5"H x 2"W x 2.5"D)
AC285-3-PM Triple
(2.25"H x 20"W x 3"D)



DESCRIPTION: Mounting Plate with single or triple broom handle holder for owner provided broom(s)



PART NUMBER: IC416-P
(1.5"W x 8"H x 1"D)

DESCRIPTION: Universal Mounting Plates with wall stand-offs for most owner provided paper towel holders with slotted mounts on back (Dispenser is not included)



PART NUMBER: IC110-P
(7.5"H x 4.5"W x 1"D)

DESCRIPTION: Universal Mounting Plate for most owner provided wide body soap / sanitizer dispensers (Dispenser and sanitizer are not included)



PART NUMBER: IC112-PM
(6"H x 5"W x 5"D)

DESCRIPTION: Mounting Plate with universal width adjustable holder for most brands of sanitizing wipes (Wipes are not included)



PART NUMBER: IC288-3-DL-PM (triple)
10"H x 18.75"W x 5"D

DESCRIPTION: Mounting Plate with triple gloves box holder, heavy duty acrylic, Dual side loading. (Glove boxes are not included)



Soap / Sanitizer / Wipes

Rail Mounting Solutions

For Every Brand and Model Number





MADE IN THE USA

PART NUMBER: IC110-P
(7.5"H x 4.5"W x 1"D)

DESCRIPTION: Universal Mounting Plate for most owner provided wide body soap / sanitizer dispensers
(Dispenser and sanitizer are not included)


NOTES: Adjustable for variable vertical placement of device



PART NUMBER: IC111-P
(7"H x 3"W x 1"D)

DESCRIPTION: Universal Mounting Plate for most owner provided narrow body soap / sanitizer dispensers
(Dispenser and sanitizer are not included)

NOTES: Dispenser model needs to be confirmed at the time of order



PART NUMBER: IC136-P
(4"H x 2.75"W x 1"D)

DESCRIPTION: Mounting Plate for owner provided Steris Kindest Care Foam Hand Sanitizer
(Dispenser and sanitizer are not included)



PART NUMBER: IC137-P
(7"H x 4"W x 1"D)

DESCRIPTION: Mounting Plate for owner provided Bobrick sanitizer dispenser
(Dispenser and sanitizer are not included)





MADE IN THE USA

PART NUMBER: IC117-P
(7"H x 2.5"W x 1"D)

DESCRIPTION: Universal Mounting Plate for owner provided aerosol foam dispensers with drip pan extension 9 oz. or 17 oz. (Steris/ Ecolab)
(Foam, dispenser and drip pan extension are not included)



PART NUMBER: IC112-PM
(6"H x 5"W x 5"D)

DESCRIPTION: Mounting Plate with universal width adjustable holder for most brands of sanitizing wipes
(Wipes are not included)



NOTES: Also works to hold some industrial style spray bottles

PART NUMBER: IC106-P
(4.5"H x 5"W x 1"D)

DESCRIPTION: Mounting Plate for owner provided Metrex Caviwipes and PDI Sani Cloth Holder
(Holder and dispenser are not included)



PART NUMBER: IC141-PM
(6.4"H x 3.3"W x 2.9" / 3.25"D)

DESCRIPTION: Mounting Plate and Holder for owner provided round bottle type dispensers / hand pumps 2.5" and 3.25" diameter



NOTES: e.g. 3M Avagard Nail Cleaner #9204
Aero Wash Eye Rinse #AW160ZEPUS,
Cederroth #725200A, Pure Nest Dispenser



MADE IN THE USA


PART NUMBER: IC122-P
(3.75"H x 3"W x 1"D)

DESCRIPTION: Mounting Plate with wall stand off
for owner provided 3M Avagard Wall
Bracket for use with 9200 / 9200C



PART NUMBER: IC145-A
(2.5"H x 5"W x 1"D)

DESCRIPTION: Adapter for owner provided
Purell 450ml manual dispenser



Paper Towel Dispensers

Rail Mounting Solutions For Every Brand and Model Number





MADE IN THE USA

PART NUMBER: IC401-A-H (horizontal)
IC401-A-V (vertical)
(2.5H x 1.5"W x 1"D)

DESCRIPTION: Universal Mounting Plates
with wall stand offs for owner
provided half size tri-fold
paper towel dispensers
(i.e. Kimberly Clark, Georgia
Pacific)
(Dispenser is not included)



PART NUMBER: IC416-P
(1.5"W x 8"H x 1"D)

DESCRIPTION: Universal Mounting Plates with wall
stand-offs for most owner provided
paper towel holders, Sharps
containers etc. with slotted mounts on
back (Dispenser is not included)

NOTES: Adapters can be moved for adjustable
vertical placement



PART NUMBER: IC405-P
(10"W x 9.125"H x 1"D)

DESCRIPTION: Mounting Plate for owner provided
GA Pacific EnMotion paper towel
dispenser (model 59462)
(Dispenser is not included)

NOTES: Specific model number must be
verified prior to order
placement





MADE IN THE USA

PART NUMBER: IC408-P
(4.5"W x 6.5"H x 1"D)

DESCRIPTION: Mounting Plate for owner provided GA Pacific SofPull paper towel dispenser (models 52109, 53009, 58201, 58217 and 58008) (Dispenser not included)

NOTES: Specific model number must be verified prior to order placement



PART NUMBER: IC411-P
(7.25"W x 1"D x 6.5"H)

DESCRIPTION: Mounting Plate for owner provided GA Pacific SofPull paper towel dispensers (models 58487 and 58470) (Dispenser not included)

NOTES: Specific model number must be verified prior to order placement



PART NUMBER: IC413-P
(7.5"W x 1"D x 7"H)

DESCRIPTION: Mounting Plate for owner provided Kimberly Clark Professional (model 09335) (Dispenser not included)

NOTES: Specific model number must be verified prior to order placement



Gloves Box Dispensers

Rail Mounting Solutions

For Every Brand and Model Number
Every Type, Size and Feed Style
Plastic, Acrylic, Metal or Wire

(STAND-ALONE MOUNTING KITS FOR OWNER SUPPLIED DISPENSERS ARE AVAILABLE UPON REQUEST)






MADE IN THE USA

PART NUMBER: IC227-1-H-PM (single)
IC227-2-H-PM (double)
IC227-3-H-PM (triple)
IC227-4-H-PM (quadruple)
Horizontal top feed application
10"H x 18.75"W x 5"D (triple)

DESCRIPTION: Mounting Plate with gloves box holder, heavy duty ABS plastic with internal dividers
(Glove boxes are not included)



PART NUMBER: IC226-1-H-PM (single)
IC226-2-H-PM (double)
IC226-3-H-PM (triple)
IC226-4-H-PM (quadruple)
Horizontal top feed application
10.125"H x 17"W x 4.25"D (triple)

DESCRIPTION: Mounting Plate with gloves box holder, ABS plastic with internal dividers.
(Glove boxes are not included)



PART NUMBER: IC288-1-DL-PM (single)
IC288-2-DL-PM (double)
IC288-3-DL-PM (triple)
10"H x 18.75"W x 5"D (triple)

IC288-4-DL-PM (Quadruple)
Vertical - side feed application only

DESCRIPTION: Mounting Plate with gloves box holder, heavy duty acrylic, Dual side loading.
(Glove boxes are not included)





MADE IN THE USA

PART NUMBER: IC297-4-PM
11.312"H x 19.75"W x 4.12"D

DESCRIPTION: Mounting Plate
with 4 square gloves box
dispenser, heavy duty acrylic,
Dual side loading.



PART NUMBER: IC287-1-H-PM (single)
IC287-2-H-PM (double)
IC287-3-H-PM (triple)
IC287-4-H-PM (quadruple)
Horizontal top feed application

IC287-1-V-PM (single)
IC287-2-V-PM (double)
IC287-3-V-PM (triple)
IC287-4-V-PM (quadruple)
Vertical side feed application

DESCRIPTION: Mounting Plate with heavy duty
acrylic gloves box holder
(Glove boxes are not included)
10"W x 18.75"H x 5"D (triple)



Please note:

Sometimes a room layout changes...

**By switching the adapters on the back,
you can easily transform your horizontal glove box
holder into a vertical glove box holder. You just need an
Allen wrench!**

Works for most Paladin glove box holders!



O.R. / Scrub Sink Solutions / PPE Dispensers

Custom Solutions

**For more Scrub Sink Solutions
please contact your local Representative!**





MADE IN THE USA

PART NUMBER: IC245-PM
(5"H x 8.75"W x 4.375"D)

DESCRIPTION: Mounting Plate with bulk dispenser,
with lid, ABS plastic and
clear PETG plastic

NOTES: Holds one box of ear loop face masks or
booties in bulk quantity



PART NUMBER: IC241-PM
(12.25"W x 11.125"H x 5.75"D)

DESCRIPTION: Mounting Plate with bulk
dispenser, with lid, double bin,
heavy duty clear acrylic



PART NUMBER: IC275-HV-PM
(5.75"H x 13.75"W x 7"D)

DESCRIPTION: Mounting Plate with PPE dispenser,
small, ABS plastic

NOTES: Universal small box dispenser for shoe covers, caps
and other items Horizontal top feed or vertical side feed application



PART NUMBER: IC276-PM
(10.125"H x 12.55"W x 10"D)

DESCRIPTION: Mounting Plate with PPE dispenser,
large, ABS plastic

NOTES: Universal large box dispenser for shoe
covers, caps and other items





MADE IN THE USA

PART NUMBER: IC271-PM
(11.625"H x 18.5"W x 4.5"D)

DESCRIPTION: Mounting Plate with PPE
dispensing cabinet, ABS plastic

NOTES: Holds two boxes of gloves, one box of face masks, and one hand sanitizer bottle



PART NUMBER: IC274-PM
(23.5"H x 17.625"W x 5"D)

DESCRIPTION: Mounting Plate with PPE dispensing cabinet,
ABS and PETG plastic

NOTES: Holds three boxes of gloves, two boxes of face masks, multi-fold towels, and patient bibs. Tabs on both sides of unit to hang patient bib clips



PART NUMBER: IC268-PM
(11.125"H x 12.25"W x 5.75"D)

DESCRIPTION: Mounting Plate with bulk dispenser,
stainless steel

NOTES: Single compartment with hinged lid holds a variety of personal protection apparel; including gloves, bouffant caps, shoe covers, etc.



PART NUMBER: IC301-PM
(15.5"W x 9"H x 4.25"D)

DESCRIPTION: Mounting Plate with dispenser,
for disposable barrier sleeves,
NA acrylic, heavy duty, clear

NOTES: Holds one box of tray sleeves, headrest covers, X-Ray sleeves or chair sleeves





MADE IN THE USA

PART NUMBER: IC542-PM

(15"H x 9.25"W x 6.25"D)

DESCRIPTION: Mounting Plate with pamphlet display, heavy duty acrylic

NOTES: With wall stand offs. Bottom of dispenser is slotted for easy cleaning



PART NUMBER: IC501-PM

(9.25"H x 7.625"D x 4"W)

DESCRIPTION: Mounting Plate with dispenser for surgical masks, heavy duty acrylic, (fits Halyard Soft Touch II, #47500) (Face masks are not included)

NOTES: With wall stand-off. Bottom of dispenser is slotted for easy cleaning



PART NUMBER: IC229-PM

(8.375"H x 5.375"W x 8.25"D)

DESCRIPTION: Mounting Plate with dispenser for bouffant cap or shoe cover, heavy duty acrylic, clear, Bowman BB-034 (Bouffant caps / shoe covers are not included)

NOTES: With wall stand off. Bottom of dispenser is slotted for easy cleaning



PART NUMBER: IC515-4-PM (4 Bin)

IC515-5-PM (5 Bin)

(7"H x 5.5"D x 23.75"W - 4 bin dispenser)

DESCRIPTION: Mounting Plate with surgical PPE dispenser, top feed, heavy duty acrylic, clear (Masks are not included)

NOTES: With wall stand off. Bottom of dispenser is slotted for easy cleaning





MADE IN THE USA

PART NUMBER: IC298-PM
(5.25"W x 5.25"D x 11"H)

DESCRIPTION: Mounting Plate with scrub brush dispenser
(Scrub brushes are not included)

NOTES: With wall stand off, heavy duty acrylic, bottom slotted for easy cleaning



PART NUMBER: IC533-PM
(Inside dimensions 5.5"H x 4.625"W x 8.125"D)

DESCRIPTION: Mounting Plate with 3M Aura surgical mask dispenser for 3M Box #1860
(Face masks are not included)

NOTES: Heavy duty acrylic, bottom slotted for easy cleaning



PART NUMBER: IC530-PM
(Inside dimensions:
8.5"H x 9 5/8"W x 3.5"D)

DESCRIPTION: Mounting Plate with 3M Aura surgical mask dispenser for 3M Box #1870
(Face masks are not included)

NOTES: With wall stand off, heavy duty acrylic, bottom slotted for easy cleaning, heavy duty acrylic



PART NUMBER: IC526-PM
(11"W x 11"H x 4"D)

DESCRIPTION: Mounting Plate with dispenser for bouffant caps for Halyard #69803
(Bouffant caps are not included)

NOTES: With wall stand-off, bottom of dispenser is slotted for easy cleaning, heavy duty acrylic





MADE IN THE USA

PART NUMBER: IC-PFM100-PM

DESCRIPTION: Mounting Plate with universal heavy duty acrylic holder for face mask with ear loops
Dispenser (Face masks are not included)

NOTES: 8.125"W x 5.625"D x 4.75"H
(dispenser only)



**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1800, 1816, 1820, 1826, 1832, 1838, 1860S

Cardinal Health: AT70021, AT71021, AT771141, AT771211

Halyard: 25867, 25868, 25869, 28797, 28802, 28806, 28808, 28810, 32856, 37525, 37525, 39117, 41802, 41805, 47080, 47085, 47090, 47107, 47108, 47117, 47119, 47127, 47295, 47297, 48100, 48105, 48201, 48215, 48296, 48310, 49214, 49310, 49700, 49701, 59927, 59928, 59929, 59931, 59936, 59938, 59939, 59940, 62113, 62115, 62356, 62357, 62362, 62363, 62364, 62367, 62369, V1025,

McKesson: 91-1000, 91-1200, 91-1400, 91-2002

MediChoice, Primagard/PriMed: 1773ELM001, 41001, 41073, 41092, 41200, 41201, 41241, 41331, 42000, 42001, 42092, 42121, 42241, 42320, 42331, 42332, 42421, 44021, 44051, 44241, 44331, 46001

Medline: NON27110, NON27120, NON27120-2, NON27122, NON27355, NON27365, NON27374, NON27375, NON27376, NON27377, NON27378, NON27385, NON27386, NON27400, NON27402, NON27408EL, NON27412, NON27412EL, NON27712EL, NON27812EL, NONSH400

Precept Prestige Ameritech/ProGear: 15100, 15101, 15130, 15150, 15155, 15160, 15200, 15211, 15212, 15230, 15240, 15260, 15300, 15320, 15420, 15500, 65-3320, 65-3322, FM81000, FM85011

Tidi: 9010, 9020, 9030, 9040



MADE IN THE USA

PART NUMBER: IC-PFM200-PM

DESCRIPTION: Mounting Plate with universal heavy duty acrylic holder for face mask with face shield dispenser
(Face masks are not included)



NOTES: 9.625"W x 4.625"D x 8.5"H
(dispenser only)

**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1120/37062, 1832, 1838R, 1870, 1870-N95, 8214, 8514, 9210/37021

Cardinal Health: AT7505-P, AT7511-WE, AT51035, AT52005, AT54635, AT71035, AT71039, AT72835, AT74535, AT74631, AT74635, AT771145, N95-ML, 22EL77

Halyard: 28804, 28809, 28821, 37505, 37510, 41803, 47095, 47137, 47147, 47148, 47149, 47298, 47500, 48225, 48227, 48228, 48237, 48247, 48248, 48257, 48298, 49210, 49212, 49230, 59926, 59930, 59937, 62114, 62116, 62368, 69801

McKesson: 91-1600, 91-2600

MediChoice, Primagard/PriMed: 41573, 41592, 42073, 42333, 42560, 42562, 43092, 43573, 44073, 45033, 45073, 46003, 46403, 77701, 77702, 77703, PG43573

Medline: NON27373A, NON27410, NON27410EL, NON27411, NON27411AG, NON27600, NON27600AG, NON27710, NON27710EL, NON27810, NON27810EL, NON27361A, NON27382, NON27408, NON27712

Precept /Prestige Ameritech/ProGear: 65-3344, 15310, 15330, 15710

Sperian: 091007, 14110452



MADE IN THE USA

PART NUMBER: IC-PFM300-PM

DESCRIPTION: Mounting Plate with universal heavy duty acrylic holder for cone style face mask dispenser
(Face masks are not included)

NOTES: 6.875"W x 5.625"D x 7.25"H
(dispenser only)



**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1800+NL, 1800NL, 1860, 1860S, 8210, 8240, 8246, 8247, 8511, 8576, 8577

Cardinal Health/ Condor: 22EL78, 22EL79, 22EL80, 22EL81, 4KMX9

Gerson: 061501-CS, 081730, 081735, 081740, 081745, 081750, 081760, 081940

Halyard: 00152, 39123, 46727, 46767, 46827, 46828, 46867, 47200, 47407, 47700, 47707, 47717, 47817, 62126, 62127, 62354, 62355, 63037, 63149, 63152

MediChoice, Primagard/PriMed: 4303, 40220

Medline: N7120

Moldex: 1200N95, 1201N95, 1207N95, 1510, 1511, 1512, 1513, 1517, 1712, 2112, 2200N95, 2200N95HV, 2201N95, 2207N95, 2210GN95, 2211GN95, 2212GN95, 2217GN95, 2250, 2300N95, 2301N95, 2307N95, 2350, 2600N95, 2601N95, 2607N95, 2700N95, 2701N95, 2707N95, 2730N100, 2731N100, 2740R95, 2741R95, 2800N95, 2801N95, 2840R95, 2841R95, 2850, 2940R95, 2941R95, 4200, 4201, 4300P95, 4600, 4601, 4700N100, 4800, D2740R95, M2600N95, M2601N95, M2700N95, M2701N95, M2707N95, M2740R95, M2741R95, M2800N95, M2801N95, M2840R95

Precept, Prestige Ameritech/ProGear: 65-3333, 65-3395, 65-3395S

Sperian: 14110388



MADE IN THE USA

PART NUMBER: IC-PFP021-PM

DESCRIPTION: Mounting Plate with
universal heavy duty
face mask dispenser holder
(Face masks are not included)

NOTES: 8.125"W x 5.625"D x 7.5"H
(dispenser only)



**FITS THE
FOLLOWING
DISPENSERS:**

Cardinal Health: AT7511, AT71021, AT771141, AT771211, AT771212,
N95-ML

Halyard: 25867, 25868, 25869, 28797, 28806, 28808, 28810, 32856, 37525,
41802, 41805, 47080, 47085, 47090, 47107, 47108, 47117, 47119, 47127,
47295, 47297, 48100, 48105, 48201, 48215, 48296, 49214, 49310, 49700,
49701, 59927, 59928, 59929, 59931, 59936, 59938, 59939, 59940, 62113,
62115, 62356, 62357, 62362, 62363, 62364

McKesson: 91-1200, 91-2002

MediChoice, Primagard/PriMed: PG41001, PG41092, PG41200, PG41201,
PG41241, PG41331, PG42092, PG42121, PG42332, 205ELM001

Medline: NONSH400, NON27110, NON27120, NON27120-2, NON27122,
NON27355, NON27365, NON27374, NON27376, NON27377, NON27378,
NON27385, NON27402, NON27408EL, NON27412, NON27412EL,
NON27600, NON27600AG, NON27710, NON27710EL, NON27712EL,
NON27810, NON27810EL, NON27812EL,

Precept: 65-3320, 65-3322, 65-3344, 15100, 15101, 15130, 15150, 15155,
15200, 15211, 15212, 15230, 15240, 15260, 15300, 15310, 15320, 15330,
15420, 15500, 15710

Prestige Ameritech/ProGear: FM81000, FM82120, FM84210, FM85011

Sperian: 14110452

Tidi: 9010, 9020, 9030, 9040



MADE IN THE USA

PART NUMBER: IC-PFP022-PM

DESCRIPTION: Mounting Plate with universal heavy duty acrylic face mask holder for face mask with ear loops dispenser
(Face masks are not included)



NOTES: 8.125"W x 5.625"D x 7.5"H
(dispenser only)

**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1832, 1838, 1900, 8214, 8514

Cardinal Health: AT7505-P, AT51035, AT71035, AT72835, AT74535, AT752005, AT771145, N95-ML

Halyard: 28800, 28804, 28809, 28820, 28821, 37505, 37510, 41803, 47137, 47147, 47148, 47298, 47500, 47625, 47650, 48207, 48208, 48210, 48220, 48225, 48227, 48228, 48237, 48247, 48248, 48257, 48297, 48298, 48390, 49210, 49212, 49215, 49216, 49230, 49235, 59926, 59930, 59937, 62114, 62116, 62368

McKesson: 91-1600, 91-2600

MediChoice, Primagard/PriMed: PG42073, PG42092, PG42121, PG42332, PG42333, PG43573

Medline: NON27361A, NON27371A, NON27378A, NON27379A, NON27382, NON27382FR, NON27408, NON27420, NON27420EL, NON27712

Precept: 15212, 15216, 15220, 15310, 15330, 15710, 65-3344

Prestige Ameritech/ProGear: FM81110



MADE IN THE USA

PART NUMBER: IC-PFP120-PM

DESCRIPTION: Mounting Plate with universal heavy duty acrylic face mask holder
for face shield dispenser
(Face masks are not included)



NOTES: 9.125"W x 4.25"D x 7.5"H
(dispenser only)

**FITS THE
FOLLOWING
DISPENSERS:**

Halyard:

00146, 28804, 28809, 28821, 41803, 47137, 47147, 47148, 47149,
47298, 47500, 48228, 48237, 48247, 48248, 48257, 48298,
59930, 59937, 62114, 62116, 69801

McKesson:

91-1600, 91-2600

MediChoice, Primagard/PriMed:

PG41573, PG41592, PG45033, PG45073

Medline:

NON27600, NON27600AG, NON27710, NON27710EL, NON27810,
NON27810EL

Precept:

65-3344, 15310, 15330, 15710

Prestige Ameritech/ProGear:

FM82120, FM84210

Sperian:

14110452

Sharps Containers

Rail Mounting Solutions For Every Brand and Model Number





MADE IN THE USA


PART NUMBER: IC309-P (locking)
(4"W x 1"D x 3.375"H)

DESCRIPTION: Mounting Plate for owner provided Daniels "Sharpsmart" sharps container (S14 or S22). (Container and Daniel provided mounting brackets are not included)




PART NUMBER: IC304-P (non-locking)
(4"W x 1"D x 3.375"H)

DESCRIPTION: Mounting Plate for owner provided Daniels "Sharpsmart" sharps container (S14 or S22). (Container and Daniel provided mounting brackets are not included)



PART NUMBER: IC303-P
(5"W x 5"H x 1"D)

DESCRIPTION: Universal Mounting Plate for owner provided Standard Kendall / Covidien Sharps Container, all sizes (Container is not included)



PART NUMBER: IC317-P
(10"W x 11"H x 1"D)

DESCRIPTION: Universal Mounting Plate for owner provided Stericycle sharps containers, all sizes (Container is not included)



Baskets and Bins

Rail Mounting Solutions

Factory Attached to Accessories – No assembly required





MADE IN THE USA

PART NUMBER: SM101-PM

(7"H x 6.75"W x 4"D)

DESCRIPTION: Mounting Plate with BP cuff basket
(Content of basket for illustration purposes only. Not included.)

NOTES: Plastic coated wire U-shaped basket. Max load 4lbs.



PART NUMBER: SM104-PM

(10"H x 16.25"W x 15.5"D /
23.5"D fully extended)

DESCRIPTION: Mounting Plate with medium deep wire bin on pull-out glide
(Content of basket for illustration purposes only. Not included.)

NOTES: Plastic coated wire basket



PART NUMBER: SM105-PM

(10"H x 23"W x 15.5"D /
23.5"D fully extended)

DESCRIPTION: Mounting Plate with large deep wire bin on pull-out glide
(Content of basket for illustration purposes only. Not included.)

NOTES: Plastic coated wire basket






MADE IN THE USA

PART NUMBER: SM109-PM
(4"H x 22"W x 5"D)

DESCRIPTION: Mounting Plate with fixed long wire basket

NOTES: Plastic coated wire basket. Max load 10lbs.



PART NUMBER: SM102-ACRYLIC
(4.25"H x 10"W x 6.5"D)
SM102-PM-WIRE
(4.25"H x 10"W x 7"D)
SM102-PM-LONG-WIRE
(4.25"H x 7"W x 10"D)
SM102-PM-HILLROM-HORIZON
(4.25"H x 10"W x 7"D)

DESCRIPTION: Mounting Plate with small basket
(Content of basket for illustration purposes only. Not included.)

NOTES: Plastic coated wire basket. Max load capacity 5lbs.



PART NUMBER: SM103-PM
(4"H x 18"W x 10"D)

DESCRIPTION: Mounting Plate with large wire basket
(Content of basket not included.)

NOTES: Plastic coated wire basket. Max load capacity 5lbs.





MADE IN THE USA

PART NUMBER: SM201-1-PM
(14"H x 12"W x 14"D)

DESCRIPTION: Mounting Plate with
1 bin tilt-out storage container
(Content of container for illustration
purposes only. Not included.)

NOTES: White plastic box with clear
plastic bins



PART NUMBER: SM201-2-PM
(14"H x 24"W x 14"D)

DESCRIPTION: Mounting Plate with
2 bin tilt-out storage container
(Content of container for
illustration purposes only. Not
included.)

NOTES: White plastic box
with clear plastic bins



PART NUMBER: SM201-3-PM
(9.5"H x 24"W x 9.5"D)

DESCRIPTION: Mounting Plate with
3 bin tilt-out storage container
(Content of container for
illustration purposes only. Not
included.)

NOTES: White plastic box with clear
plastic bins





MADE IN THE USA

PART NUMBER: SM201-4-PM
(8.25"H x 23.75"W x 5.5"D)

DESCRIPTION: Mounting Plate with
4 bin tilt-out storage container
(Content of container for
illustration purposes only. Not
included.)

NOTES: White plastic box with clear
plastic bins



PART NUMBER: SM201-5-PM
(6.5"H x 23.75"W x 4.5"D)

DESCRIPTION: Mounting Plate with
5 bin tilt-out storage container
(Content of container for
illustration purposes only. Not
included.)

NOTES: White plastic box with clear
plastic bins



PART NUMBER: SM291
(12"W x 11"H x 0.5"D)

DESCRIPTION: Mounting Plate for owner provided
Teleflex Arrow EZ-iO cabinet.



Clips, Hooks and Wire Management

Rail Mounting Solutions





MADE IN THE USA

PART NUMBER: SM301-PM-1 (1" hook)
SM301-PM-2 (2" hook)
SM301-PM-4 (4" hook)
(2.5"H x 1.5"W x X"D
length of hook ¾", 1 ½", 2" and 4")

DESCRIPTION: Mounting Plate with 1 hook hanger
powder coated metal wire hook

NOTES: Hook and integrated mounting plate in a powder coat painted anti-microbial finish(AMF)



PART NUMBER: SM302-PM-1 (1" hook)
SM302-PM-2 (2" hook)
SM302-PM-4 (4" hook)
(2.5"H x 12"W x X"D
length of hooks available: ¾", 2", and 4")

DESCRIPTION: Mounting Plate with 3 hook hanger
powder coated metal wire hooks

NOTES: Hooks and integrated mounting plate in a powder coat painted anti-microbial finish (AMF)



PART NUMBER: SM304-PM
(3.5"H x 1.5"W x 0.75"D)

DESCRIPTION: Adapter with cord bundler,
ABS antimicrobial

NOTES: Slide action loop and clip to secure cords, cables, and wires, or to hold stethoscope.





MADE IN THE USA


PART NUMBER: SM303-PM
(1.5"H x 1.5"W x 1.75"D)

DESCRIPTION: Mounting Plate
with lead apron hook
for owner provided lead
aprons



PART NUMBER: SM308-PM
(1.5"H x 1.5"W x 1.75"D)

DESCRIPTION: Mounting Plate with
stationary double hook,
stainless steel, for owner
provided hangers



PART NUMBER: SM307-PM
(1.75"H x 3.25"W x 1"D)

DESCRIPTION: Mounting Plate with
swing / adjustable triple hook,
stainless steel



PART NUMBER: SM309-PM
(1.5"H x 1.5"W x 1.75"D)

DESCRIPTION: Mounting Plate
with U-shape universal Hook,
Heavy duty
Can be used with owner provided
Kendell SCD700 series device





MADE IN THE USA

PART NUMBER: SM401-PM-3
(2.5"H x 9.75"W x 1"D)

DESCRIPTION: Mounting Plate
with 3 clip hangers
(Items on hangers for illustration purposes only. Not included.)

NOTES: Metal clips on anti-microbial (AMF) powder coat painted back-plate



PART NUMBER: SM401-PM-6
(2.5"H x 22"W x 1"D)

DESCRIPTION: Mounting Plate with 6 clip hangers
(Items on hangers for illustration purposes only. Not included.)

NOTES: Metal clips on anti-microbial (AMF) powder coat painted back-plate



PART NUMBER: SM406-PM
(1.5"H x 13"W x 1"D)

DESCRIPTION: Mounting Plate
with document holder
(Items on holder for illustration purposes only. Not included.)


NOTES: Friction fit paper display rail.
Release document by lifting bottom of sheet up towards ceiling.



Shelves and Work Surfaces

Rail Mounting Solutions

Factory Attached to Accessories – No assembly required





MADE IN THE USA

- PART NUMBER:** **S201-PM**
 (10"W x 8"D)
 S203-PM
 (14.5"W x 13.5"D)
- DESCRIPTION:** **Mounting Plate with fixed shelf**
 with lip edges
 (Monitor not included)
- NOTES:** Metal with anti-microbial (AMF)
 powder coat painted finish.
 Maximum load 5lbs.
- Includes 0.75" lip on 3 sides



- PART NUMBER:** **S205-PM**
 (31"W X 10"D)
- DESCRIPTION:** **Mounting Plate**
 with vented shelf, stainless steel
 Maximum load 5lbs



- PART NUMBER:** **S242-PM small (15"W x 10"D)**
 S241-PM medium (20"W x 13"D)
 S244-PM large (24"W x 18"D)
- DESCRIPTION:** **Mounting Plates**
 with heavy duty flip shelf
 (Item on shelf for illustration
 purposes only. Not included.)
- NOTES:** Scratch resistant food grade high
 density polyethylene.
 Indented spill lip.
 Up to 100lbs load capacity.





MADE IN THE USA

PART NUMBER: CF201-P

DESCRIPTION: Mounting Plate
for Herman Miller C-Frame

NOTES: For quick and stable mounting
of C-Frame to the “Evolution”
Fairfield-style rail and to the
Paladin “M.T.S” Mobile
Treatment Station

C-Frame unit sold separately

C-Frame unit customizable



PART NUMBER: CF101-K

DESCRIPTION: Mounting Plate with side rail kit
for Herman Miller C-Frame
(Items on cart for illustration purposes
only. Not included.)

NOTES: Add equipment rail to Herman
Miller C-Frame for easy
mounting of accessories to side
of C-Frame unit.

2 rails included.

C-Frame unit sold separately



Monitor Mounting and Cable Management Tools

Smart Solutions Workspace Flexibility





MADE IN THE USA

PART NUMBER: GCX101-19-2 (19"H x 4"W x 2"D)
GCX101-25-2 (25"H x 4"W x 2"D)
GCX101-31-2 (31"H x 4"W x 2"D)
GCX101-21-2 IntegrisHW

**Additional custom lengths
available**

DESCRIPTION: GCX101-19-2 channel to span (2) rails
at 7-16" between rails
GCX101-25-2 channel to span (2) rails
at 12-22.5" between rails
GXC101-31-2 channel to span (2) rails
at 15-28" between rails
GCX101-21-2-IntegrisHW
(21"H x 2.85"W x 3"D)

(Monitor and BP cuff not included)

NOTES: To mount a variety of equipment, arms, and management tools. Includes (2) pre-installed adapters.



PART NUMBER: GCX101-49-3 (49"H x 4"W x 2"D)
Additional custom lengths available

DESCRIPTION: GCX channel to span (3) rails
at 10-20" between rails, but no more than 35"
maximum and 25" minimum total between top
and bottom rail. (Monitor, swingarm and
keyboard are not included)

NOTES: To mount a variety of equipment such as
computer arms and equipment or swing-arm
IV pole. Includes (3) pre-installed adapters.





MADE IN THE USA

PART NUMBER: GCX149 (5.5"H x 5.5"W x 2"D)

DESCRIPTION: GCX channel for owner provided GCX style swing arms (Monitor and swing arm (PR-GCX-FLP-0009-02) are not included)


NOTES: Integrates the GCX channel to a single equipment rail. GCX channel is factory attached to rail mounting adapter.



PART NUMBER: DEM163-PM
(6"H x 3.5"W x 8"D)

DESCRIPTION: Mounting Plate with articulating / tilting VESA mount arm (Monitor is not included)

NOTES: 75/100mm with rail mount kit



PART NUMBER: AC126-PM-2 (double)
(8.5"H x 5.5"W x 2.5"D)
AC126-PM-3 (triple)
(8.5"H x 9"W x 2.5"D)
AC126-PM-4 (quadruple)
(8.5"H x 12"W x 2.5"D)
AC126-PM-5 (quintuple)
(8.5"H x 15"W x 2.5"D)
AC126-PM-6 (sextuple)
(8.5"H x 18"W x 2.5"D)

DESCRIPTION: Mounting Plate with Cord Caddy - cord, hose storage and retrieval system
(All cables on the cord caddy are for illustration purposes only and not included)





MADE IN THE USA

PART NUMBER: AC109-PM
(6"H x 1.5"D x 3"W)

DESCRIPTION: Mounting Plate with the "Hold Everything" mount for owner provided tablets and smartphones...

NOTES: User adjustable width to accommodate most mobile devices, computer accessories, hand-held controllers
Includes one pair of brackets (left and right)



PART NUMBER: AC125-PM
(3"W x 2.5"H x 1.5"D)

DESCRIPTION: Mounting Plate with holder for owner provided pillow speaker/nurse call (Crest Healthcare and Curbell Medical Style) / patient angled phone handset (Med-Pat)

NOTES: Shown with SM304 cord bundler (sold separately)



PART NUMBER: AC185-PM
(3.25"H x 3"W x 2"D)

DESCRIPTION: Mounting Plate with nurse call cord button holder for owner provided nurse call

NOTES: Shown with SM304 cord bundler (sold separately)





MADE IN THE USA

PART NUMBER: AC285-1-PM Single
(2.5"H x 2"W x 2.5"D)

DESCRIPTION: Mounting Plate with single broom handle holder
for owner provided broom



PART NUMBER: AC285-3-PM Triple
(2.25"H x 20"W x 3"D)

DESCRIPTION: Mounting Plate with triple broom handle holder
for owner provided brooms



PART NUMBER: AC121-A
(2.5"H x 1"D x 1.5"W)

DESCRIPTION: Adapter for owner provided
Posey Chair Alarm, Fall
Monitor or Sitter Bracket



PART NUMBER: AC233-A
(2 3/4 "W x 4"H x 1"D)

DESCRIPTION: D2HT Metal Adapter
for owner provided Stanley
M200 Chair Fall Alarm





MADE IN THE USA

PART NUMBER: SM304-PM
(3.5"H x 1.5"W x 0.75"D)

DESCRIPTION: Adapter with cord bundler,
ABS antimicrobial
(Stethoscope for illustration
purposes only. Not included.)

NOTES: Slide action loop and clip to secure
cords, cables, and wires, or
to hold stethoscope.



PART NUMBER: SM311-B-PM backward facing
SM311-F-PM forward facing
(1.5"W x 1.25"D x 3.5"H)

DESCRIPTION: Mounting Plate with cable/gas
hose management adapter
for owner provided hoses


NOTES: Accommodates cables and gas
hoses for secure and orderly
management
Powder coat anti-microbial
finish (AMF)



PART NUMBER: SM306-PM
(8.25"W x 1"D x 1.5"H)

DESCRIPTION: Mounting Plate with cord cleat
(Cable shown for illustration
purposes only. Not included.)

NOTES: Easy management of cables and
cords



Medical Gas and Vacuum Management





Vacuum Bottle Slides

Do not include vacuum bottles



PART NUMBER: AC101-PM-FAIR
(2.125"W x 1"D x 4.125"H)

DESCRIPTION: Mounting Plate with vacuum bottle slide for Fairfield style rail, with anti microbial adapter. Will also fit Hillrom Integris Headwall System.



PART NUMBER: AC101-PM-HILLROM-HORIZON
(2.125"W x 1"D x 4.125"H)

DESCRIPTION: Mounting Plate with vacuum bottle slide for HillRom Horizon Style Headwall, horizontal, with anti microbial adapter



PART NUMBER: AC101-PM-HILLROM-VERT
(2.125"W x 1"D x 4.125"H)

DESCRIPTION: Mounting Plate with vacuum bottle slide for HillRom Vertical Channel Mount, with anti microbial adapter



PART NUMBER: AC101-PM-MODULAR
(2.125"W x 1"D x 4.125"H)

DESCRIPTION: Mounting Plate with vacuum bottle slide for Modular and Amico Horizontal Headwalls and Rail Systems, with anti microbial adapter



PART NUMBER: AC101-PM-SWIV
Available for Fairfield Style Rail OR HillRom Vertical Mount
Please specify at time of order!
(2.125"W x 1"D x 4.125"H)

DESCRIPTION: Mounting Plate with vacuum bottle slide and Fairfield Swivel Mount OR Hillrom Vertical Channel Swivel Rail Mount, with anti microbial adapter





MADE IN THE USA

PART NUMBER: GBA105-PM (Air)
GBO105-PM (Oxygen)
GBV105-PM (Vacuum)
(4"W x 3"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical-gas block
with 2 x 1/8" NPT male front outlets.
Air and oxygen blocks come with an
additional 90 degrees DISS center outlet
check valve. Vacuum block does not include
center check valve.

NOTES: Flow meters and vacuum regulator not included.



PART NUMBER: GBA106-PM (Air)
GBO106-PM (Oxygen)
GBV106-PM (Vacuum)
(4"W x 3"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical gas block
with 2 x DISS male front check valves. Air and
oxygen blocks come with an additional 90
degrees DISS center outlet check valve. Vacuum
block does not include center check valve.

NOTES: Flow meters and vacuum regulator not included.



PART NUMBER: GBA110-PM (Air)
GBO110-PM (Oxygen)
GBV110-PM (Vacuum)
(4"W x 3"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical gas block
with 2 x female quick connect front outlets. Available with
Ohio, Chemetron, Oxyquip and Schraeder fittings. Air and oxygen
blocks come with an additional 90 degrees DISS center outlet check
valve. Vacuum block does not include center check valve.

NOTES: Flow meters and vacuum regulator not included.






MADE IN THE USA

PART NUMBER: GBA104-PM (Air)
GBO104-PM (Oxygen)
GBV104-PM (Vacuum)
(4"W x 3"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical gas block
with 1 x 1/8" NPT male fitting front outlet,
1 x female quick connect front outlet. Available
with Ohio, Chemetron, Oxyquip and Schraeder
fittings. Air and oxygen blocks come with an
additional 90 degrees DISS center outlet check
valve. Vacuum block does not include center
check valve.

NOTES: Flow meters and vacuum regulator not included.



PART NUMBER: SBO201-PM (Oxygen)
SBA201-PM (Air)
SBV201-PM (Vacuum)
(3"W x 2"H x 1.5"D)

DESCRIPTION: Mounting Plate with wall outlet splitter block
DISS or Quick connect to 2 DISS check valves.
Available with DISS Ohio, Chemetron, Oxyquip
and Schraeder fittings.



PART NUMBER: AC104-PM
(2"W x 1"D x 3"H)

DESCRIPTION: Mounting Plate with universal
blender slide, for any owner
provided brand of blender.
Blender not included.



NOTES: Bracket slide with anti-microbial powder coat finish.



MADE IN THE USA

PART NUMBER: VM201-PM
(13"W x 3.5"H x 2.5"D)

DESCRIPTION: Mounting Plate with double ventilator manifold with DISS outlets:
2 air and 2 oxygen

NOTES: Integrated blender mount is standard. Blender is not included. Available with OHIO, Chemetron, Quick Connect or Puritan Bennet type fittings



PART NUMBER: OTH101-PM
(16"H x 5.5"W x 5.5"D)

DESCRIPTION: Mounting Plate with wire metal oxygen cylinder holder
for D or E tank cylinders

NOTES: Does not include cylinder



PART NUMBER: AC163-PM (Plate with suction pump)
AC163-P (Plate only)

DESCRIPTION: PM- Mounting Plate with LCSU 4 Laerdal suction pump, (7"W x 9"H x 8"D)

NOTES: P- Mounting Plate for owner provided LCSU 4 Laerdal suction pump,
(1.5"W x 8"H x 1"D)





MADE IN THE USA

PART NUMBER: FN508-PM (Plate with oxygen station)
FN508-P (Plate only)

DESCRIPTION: PM - Mounting Plate with FERNO #T70710-N oxygen station and 2 carbon fiber oxygen tanks (18"W x 34.5"H x 10"D)
P - Mounting Plate for owner provided #T70710-N FERNO oxygen station and tanks


NOTES: Back wall drilled for Paladin adapter mounting



PART NUMBER: GBO115-PM
(18"W x 5"H x 7.5"D)

DESCRIPTION: Mounting Plate with FERNO Multi Manager Gas Block with multiple gas flow regulators (8)


NOTES: Primarily use for mass casualty, carts and free-standing tripod



PART NUMBER: GBO116-PM
(8"H x 6.5"W x 5"D)

DESCRIPTION: Mounting Plate with FERNO Mini Manager Compact Gas Block with multiple gas flow regulators (5)

NOTES: Primarily use for mass casualty, carts and free-standing tripod



PART NUMBER: GBO114-PM
(7"W x 3"H x 2.5"D)

DESCRIPTION: Mounting Plate with FERNO Micro Manager 5-way DISS oxygen gas block

NOTES: Primarily use for mass casualty, carts and free-standing tripod





PART NUMBER: GHA201-2-90 (2ft – DISS female to DISS female)
GHA201-3-90 (3ft – DISS female to DISS female)
GHA201-4-90 (4ft – DISS female to DISS female)
GHA202-2-90 (2ft – DISS to Chemetron Style Quick Connect)
GHA202-3-90 (3ft – DISS to Chemetron Style Quick Connect)
GHA202-4-90 (4ft – DISS to Chemetron Style Quick Connect)
GHA203-2-90 (2ft – DISS to Ohio Style Male Quick Connect)
GHA203-3-90 (3ft – DISS to Ohio Style Male Quick Connect)
GHA203-4-90 (4ft – DISS to Ohio Style Male Quick Connect)

DESCRIPTION: Air hose with 90-degree elbow fitting, clear nonconductive.

PART NUMBER: GHO201-2-90 (2ft – DISS female to DISS female)
GHO201-3-90 (3ft – DISS female to DISS female)
GHO201-4-90 (4ft – DISS female to DISS female)
GHO202-2-90 (2ft – DISS to Chemetron Style Quick Connect)
GHO202-3-90 (3ft – DISS to Chemetron Style Quick Connect)
GHO202-4-90 (4ft – DISS to Chemetron Style Quick Connect)
GHO203-2-90 (2ft – DISS to Ohio Style Male Quick Connect)
GHO203-3-90 (3ft – DISS to Ohio Style Male Quick Connect)
GHO203-4-90 (4ft – DISS to Ohio Style Male Quick Connect)

DESCRIPTION: Oxygen hose with 90-degree elbow fitting, clear nonconductive.

PART NUMBER: GHV201-2-90 (2ft – DISS female to DISS female)
GHV201-3-90 (3ft – DISS female to DISS female)
GHV201-4-90 (4ft – DISS female to DISS female)
GHV202-2-90 (2ft – DISS to Chemetron Style Quick Connect)
GHV202-3-90 (3ft – DISS to Chemetron Style Quick Connect)
GHV202-4-90 (4ft – DISS to Chemetron Style Quick Connect)
GHV203-2-90 (2ft – DISS to Ohio Style Male Quick Connect)
GHV203-3-90 (3ft – DISS to Ohio Style Male Quick Connect)
GHV203-4-90 (4ft – DISS to Ohio Style Male Quick Connect)

DESCRIPTION: Vacuum hose with 90-degree elbow fitting, clear nonconductive.

Infusion Management and Lighting





MADE IN THE USA

PART NUMBER: IM301-K
(37"H x 5"D x 2"W)

DESCRIPTION: Rail Mounting Kit with IV swing-arm holder and 4-hook IV bag hanger

NOTES: Infusion pump mounting swing arm for attachment of infusion pump and other pole-mount devices
Includes adjustable telescoping IV pole with 4 bag hanger

Does not include GCX channel



PART NUMBER: IM200-K-4-3L
(36" x 20" Pole x 4"D)

DESCRIPTION: Rail mounting kit with IV pole dual rail solution, 4" offset from rail, consists of 2x rail mount IV pole holder and 1 IV pole, 3 hooks (IV bags are not included)



PART NUMBER: IM401
(1"Diameter Pole x 5"H on 3" deep post mount)

DESCRIPTION: Mounting Plate with universal post mount for any device with a c-clamp pole mount

NOTES: For Rail Attachment for Pole Mounted Equipment

*Available with custom stand-off depth (greater or fewer than 5") as defined by customer.
Available for HillRom.*





MADE IN THE USA

DESCRIPTION: **Mounting Kit for Waldmann
Triango 100 minor procedure light**

NOTES: Conversion kit to mount
Waldmann light wall mount
system to Paladin equipment rail.

Available complete with
Waldmann light for turn-key
installation.

Please call our office for details.



PART NUMBER: **LM101-PM**
 (2.5"H x 3"D x 1.5"W)
 LM102-PM (GCX channel mount kit)
 (3"H x 3"D x 2.5"W)

DESCRIPTION: **LM101-PM - Mounting Plate with
Waldmann LED flex arm light mount for
Owner provided Waldmann light
(shown)**

**LM102-PM - Mounting Plate with
Waldmann LED flex arm light mount and
GCX slide plate for owner provided
Waldmann light (shown)**

NOTES: Anti-microbial powder coat painted
finish (AMF) machined metal post
with integrated allen screw to
secure light post. Available
complete with Waldmann light for
turn-key installation please call our
office for details.



Miscellaneous Accessories



Room #:	Today's Date:
Phone #:	
Preferred Name:	
Physician:	
Care Team	
Nurse:	
Charge Nurse:	
Tech:	
Patient Education:	

PALADIN HEALTHCARE VuEase®





MADE IN THE USA

PART NUMBER: Tool Set

DESCRIPTION: Hex Driver Bit Set

NOTES: For installation and maintenance
of Paladin Healthcare rail adapters



PART NUMBER: AC174-4-PM (4 slot) HIPAA (12.5"W x 19.5"H x 3.5"D)
AC174-3-PM (3 slot) HIPAA (12.5"W x 15.75"H x 3.5"D)
AC174-2-PM (2 slot) HIPAA (12.5"W x 12.375"H x 3.5"D)
AC174-1-PM (1 slot) HIPAA (12.5"W x 8"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical chart holder, HIPAA,
horizontal

PART NUMBER: AC174-4-PM (4 slot) (12.5"W x 19.5"H x 3.5"D)
AC174-3-PM (3 slot) (12.5"W x 15.75"H x 3.5"D)
AC174-2-PM (2 slot) (12.5"W x 12.375"H x 3.5"D)
AC174-1-PM (1 slot) (12.5"W x 8"H x 3.5"D)

DESCRIPTION: Mounting Plate with medical chart holder, clear,
horizontal



PART NUMBER: AC159-PM-16x26
(16"W x 26"L x 1"D)
AC159-PM-14x18
(14"W x 18"L x 1"D)
AC159-PM-12x15
(12"W x 15"L x 1"D)

DESCRIPTION: Mounting plate with
heavy-duty acrylic single patient information
board

NOTES: Insert is not included





MADE IN THE USA

PART NUMBER: **AC173-4-PM-HIPAA (4 slot)**
(9.5"W x 21.625"H x 3.5"D)
AC173-3-PM-HIPAA (3 slot)
(9.5"W x 19.5"H x 3.5"D)
AC173-2-PM-HIPAA (2 slot)
(9.5"W x 14"H x 3.5"D)
AC173-1-PM-HIPAA (1 slot)
(9.5"W x 10"H x 3.5"D)

DESCRIPTION: **Mounting Plate with medical chart holder, vertical**

Acrylic, opaque front face, HIPAA compatible, with adapter mount



PART NUMBER: **AC173-4-PM (4 slot)**
(9.5"W x 21.625"H x 3.5"D)
AC173-3-PM (3 slot)
(9.5"W x 19.5"H x 3.5"D)
AC173-2-PM (2 slot)
(9.5"W x 14"H x 3.5"D)
AC173-1-PM (1 slot)
(9.5"W x 10"H x 3.5"D)

DESCRIPTION: **Mounting Plate with acrylic chart holder clear front / magazine rack, vertical**



PART NUMBER: **AC149-PM-14.75x12 (14.75"W x 12"L x 2"D)**
AC149-PM-16x26 (16"W x 26"L x 2"D)
AC149-PM-17.75x15 (17.75"W x 15"L x 2"D)

DESCRIPTION: **Mounting Plate with Vu Ease Patient Information Board**

NOTES: The information boards are built with high-quality tempered glass and include a limited lifetime warranty. Insert is not included.





MADE IN THE USA

PART NUMBER: AC196-PM
(16.25"W x 19"L x 2"D)

DESCRIPTION: Mounting Plate with
24 slot business card holder



PART NUMBER: IC546-PM
(5.25"W x 5.75"L x 2"D)

DESCRIPTION: Mounting Plate with
"small Kleenex" paper tissue
dispenser and wall stand off,
clear acrylic



PART NUMBER: IC512-PM
(4¾"W x 9¾"L x 2"D)

DESCRIPTION: Mounting Plate with
with "Georgia Pacific Box" tissue
dispenser and wall stand off,
clear acrylic



Notes: To fit box: 4¾"W x 9¾"H x 2"D

PART NUMBER: AC148-PM (17.75"W x 13.5"L)

DESCRIPTION: Mounting Plate with
contemporary style mirror



PART NUMBER: AC106-PM
(6"W x 7"H x 1"D - plate)

DESCRIPTION: Mounting Plate with 12 Qt. small vanity waste basket (Rubbermaid)

NOTES: Waste basket available in white or beige



PART NUMBER: IC519-PM
(8"W x 1"D x 6"H - plate)

DESCRIPTION: Mounting Plate with TIDI PPE cabinet (Content of cabinet for illustration purposes only. Not included.)



PART NUMBER: SM210-PM
(8-3/4W x 7-1/4H x 6D)

DESCRIPTION: Mounting Plate with single door narcotic box, dual lock

NOTES: Rail mount kit and double lock narcotic box, 20-gauge steel, beige



PART NUMBER: IC319-P
(1.5"W x 1"D x 4.25"H)

DESCRIPTION: Mounting Plate for Emesis bag dispenser

NOTES: Bags and dispenser are not included





MADE IN THE USA

PART NUMBER: SM220-PM
(23.5"W x 11.75"H x 11.75"D)

DESCRIPTION: Mounting Plate
with white supply cabinet with
¼" clear polycarbonate door



PART NUMBER: SM221-PM
(23.5"W x 11.75"H x 11.75"D)

DESCRIPTION: Mounting Plate
with supply cabinet with
¼" frosted polycarbonate door





MADE IN THE USA

PART NUMBER: SM206-PM

(6"W x 4"D x 18"H)

DESCRIPTION: Mounting Plate with catheter container mount, for owner provided catheters

NOTES: Includes catheter container with 2 internal dividers and pre-installed adapter mount.



PART NUMBER: A112-P

(2.5"H x 1.5"W x 1"D)

DESCRIPTION: Mounting Plate for owner provided CareFusion Razor



PART NUMBER: DEM157-PM (4 3/4"H x 5 3/8" Dia)

DESCRIPTION: Mounting Plate with antimicrobial coated cup for owner provided pens, pencils, tongue depressors, etc.

NOTES: White, powder coated



PART NUMBER: DEM130-PM-2 (2ft cord)

DEM130-PM-6 (6ft cord)

DEM130-PM-10 (10ft cord)

DEM130-PM-15 (15 ft cord)

DEM130-P (Plate only)

(9"H x 2.5"W x 1"D - plate)

DESCRIPTION: PM - Mounting Plate with Tripp-Lite hospital grade power strip, 4 outlet power strip

P- Mounting Plate for owner provided Tripp-Lite hospital grade power strip





MADE IN THE USA

PART NUMBER: FN501-PM (FERNO Safe Pac Holder)
(11"H x 13"W x 1"D)
PR-FN506B/R (medium bag)
PR-FN501 (large bag)

DESCRIPTION: FN501-PM - Mounting Plate with FERNO Safe-Pac Bag Holder for owner provided medium or large FERNO bags

NOTES: Bags available in 3 different colors: red, white and blue



PART NUMBER: FN101-PM-3
(16"W x 18"H x 2.5"D)

DESCRIPTION: Mounting Plate with FERNO Pouch Mounting Plate, incl. wall standoff for owner provided FERNO supply pouches

NOTES: For micro, small and large FERNO supply pouches



PART NUMBER: PR-FN-502-MED (medium pouch)
(4.5"H x 8"W x 2.5"D)
PR-FN-503-PM (large pouch)
(10"H x 8"W x 2.5"D)

DESCRIPTION: FERNO Inventory Supply Pouch



PART NUMBER: FN502-PM
(8"W x 4"H x 1.5"D)

DESCRIPTION: Mounting Plate with FERNO Pouch mounting plate

NOTES: This plate is for every size of FERNO Pouch



Architectural Headwall Systems



Applications:

ICU

NICU

PACU

Transitional

Step-Down

Patient Rooms

LDR Post-Partum

Key Features

- Easily configurable, scalable and adaptable
- Ergonomic placement of services to minimize caregiver bending and reaching
- Integrated workflow components to facilitate charting, eliminates unnecessary steps and simplifies bedside patient care
- Flexible storage and equipment organization options; compartments that adapt to clinician needs
- Easy to clean, reveal-less and faceplate-plate design
- Streamlined installation; headwall is assembled at the factory and shipped pre-wired and pre-plumbed with integral leveling system



Applications:

ICU

NICU

PACU

Transitional

Step-Down

Patient Rooms

LDR Post-Partum

Key Features

- Easily configurable, scalable and adaptable
- Ergonomic placement of services to minimize caregiver bending and reaching
- Integrated workflow components to facilitate charting, eliminates unnecessary steps and simplifies bedside patient care
- Flexible storage and equipment organization options; compartments that adapt to clinician needs
- Easy to clean, reveal-less and faceplate-plate design
- Streamlined installation; headwall is assembled at the factory and shipped pre-wired and pre-plumbed with integral leveling system



Applications:

ICU

NICU

PACU

Transitional

Step-Down

Patient Rooms

LDR Post-Partum

Key Features

- Easily configurable, scalable and adaptable height and width

- Ergonomic placement of services to minimize caregiver bending and reaching

- Integrated Equipment Management Rail System to facilitate change and optimizes bedside workflow

- May be recessed with only 3rd depth in the wall and less than 1st reveal on the surface. May be surface mounted as well.

- May incorporate back-lit artwork and high-style finishes.



Applications:

ICU

NICU

PACU

Transitional

Step-Down

Patient Rooms

LDR Post-Partum

Key Features

- Easily configurable, scalable and adaptable
- Ergonomic placement of services to minimize caregiver bending and reaching
- Integrated workflow components to facilitate charting, eliminates unnecessary steps and simplifies bedside patient care
- Flexible storage and equipment organization options; compartments that adapt to clinician needs
- Easy to clean, reveal-less and faceplate-plate design
- Streamlined installation; headwall is assembled at the factory and shipped pre-wired and pre-plumbed with integral leveling system



Applications:

ICU

NICU

PACU

Transitional

Step-Down

Patient Rooms

LDR Post-Partum

Key Features

- Easily configurable, scalable and adaptable

- Ergonomic placement of services to minimize caregiver bending and reaching

- Integrated workflow components to facilitate charting, eliminates unnecessary steps and simplifies bedside patient care

- Flexible storage and equipment organization options; compartments that adapt to clinician needs

- Easy to clean, reveal-less and faceplate-plate design

- Streamlined installation; headwall is assembled at the factory and shipped pre-wired and pre-plumbed with integral leveling system



Mobile Treatment and Transport Systems



PART NUMBER: MC101

DESCRIPTION: F.T.S.
Flexible Treatment Station

- NOTES:**
- Manage all clinical and support equipment in forward triage or hallway overflow conditions
 - Use as a flu screening cart in the ED waiting areas to alleviate throughput
 - Able to supply locking drawers and a variety of work surfaces and file organizers for admission forms, etc.
 - Stable weighted base that is powder-coated in an anti-microbial finish
 - Heavy-duty rubberized non-skid casters, 2 of which are lockable
 - (20) Linear Feet of usable equipment rail space for custom location of all accessories based on immediate need
 - Adjustable transport docking device will adapt to uneven surfaces or emergency Trendelenburg positioning
 - Includes (2) infusion poles and (5) EVOLUTION equipment rails
 - 26"W x 26"D x 60"H (without IV pole) – 77"H with max IV pole





PART NUMBER: SM501-4-K

DESCRIPTION: Small Cart with 4 rails including
1x SM109 (wide basket)
1x SM104 (deep slide basket)
6x SM102 (small long basket)
1x IC288-3-DL-PM (triple glove box, dual side load)
1x IC110-P (Mounting plate for owner provided sanitizers)
1x IC112-P (Mounting plate for owner provided sanitizers)
Base configuration comes with food grade top shelf as work surface.

- NOTES:**
- Manage all clinical and support equipment in forward triage or hallway overflow conditions
 - Use as a flu screening cart in the ED waiting areas to alleviate throughput
 - Able to supply locking drawers and a variety of work surfaces and file organizers for admission forms, etc.
 - Stable weighted base that is powder-coated in an anti-microbial finish
 - Adjustable transport docking device will adapt to uneven surfaces or emergency Trendelenburg positioning



PART NUMBER: MC102

DESCRIPTION: Surge / Mass Casualty Rapid Deployment
Equipment Management System

NOTES: Portable and easily stored Equipment Management System for rapid deployment in any situation. Application on rough terrain is accommodated through use of heavy duty Tri-Pod Stand which will “self level”.

- Includes (2) tri-pods with (2) carrying cases (2) 4ft rails with rail joiner in canvas pouch
- Sturdy and stable
- Pre-Hospital or In-Hospital
- Perfect for MCI Tent Applications
- Can be stored in a closet, MCI Trailer or heavy rescue type vehicle
- Delivers Med Gas, mounts PPE and monitoring devices on rail and off the floor
- Complete Integration with FERNO Intraxx components



PART NUMBER: IM601

**DESCRIPTION: METS 1:
Mobile Equipment Transportation System**

NOTES: Bed mountable equipment transport solution to facilitate patient monitoring and infection control right at the bed during patient transportation or for temporary hallway use.



PART NUMBER: IM501

**DESCRIPTION: METS 2:
Mobile Equipment Transportation System**

- NOTES:**
- Patient transport cart which docks to a patient bed or stretcher to manage and organize equipment during patient transport
 - Includes (2) infusion poles and (5) EVOLUTION equipment rails
 - Supports medical delivery, patient monitoring, supply, and infection control tools
 - Stable weighted base that is powder-coated in an anti-microbial finish
 - Heavy-duty rubberized non-skid casters, 2 of which are lockable
 - (20) Linear Feet of usable equipment rail space for custom location of all accessories based on immediate need
 - Adjustable transport docking device will adapt to uneven surfaces or emergency Trendelenburg positioning





Critical Care Anywhere Quickly Convert Alternate Care Sites



We are challenged with providing APPROPRIATE critical care environments in nontraditional spaces. Paladin Healthcare LLC can deliver total solutions for fixed and mobile conditions, with a VERY short installation time. Self-contained high-volume oxygen systems, battery powered suction devices and all aspects of PPE managed in small footprint mobile or rail mounted systems.

Visit our COVID-19 SURGE page on our web site www.PaladinHC.com for more information.





MADE IN THE USA


EQUIPMENT MANAGEMENT



CHALLENGE:

Mounting equipment directly to the wall is a costly proposition when changes need to be made.

The inevitable changes damage finishes and furnishings that often cannot be repaired—Corian, wall covering, laminate, veneer, etc.



SOLUTION:

Create and support a highly functional clinical environment now and make changes for decades to come without lost patient room revenue.

The “Evolution” Medical Equipment Management Rail System provides for nearly unlimited reconfiguration of equipment and accessories. Simply unlock the adapter and items snap-on and snap-off at will.

Easily make any necessary changes while the patient is in the room and without compromising their safety or comfort.

- Compatible with large scale equipment management tools from GCX, Human Scale, Ergotron, and many others
- Make changes without moving patients. No lost patient stay revenue.
- Complete new product changes in days instead of weeks or months.
- Proprietary embedded anti-microbial technology in our plastic adapters and powder coated items.

Optimize your Patient Care Environment

CHALLENGE:

The day after open house, medical facilities are besieged with clinical operations realities.

Often times the critical management tools were “Value Engineered” out of the project, leaving the staff to work around and “deal with” what’s in place. They have no means to define and manage the space and equipment they need.



SOLUTION:

Utilizing the “Evolution” Medical Equipment Management Rail System empowers the Clinical Staff to easily adapt and make changes to items normally fixed in place. They can create and define an optimum working environment without additional construction.

Furthermore, a neat, organized environment creates a positive impression on the patient’s and family’s perception of care.

Optimized patient care is a win-win-win for user teams, patients, and administration.

- Create user defined space now and for decades of changes to come
- Attach and move clinical tools including monitors and computer equipment without damage to the walls
- Multiply and change locations of medical gasses
- Support LEAN initiatives and manage inventory control
- Clear floor space at the bedside with IV pole management


Transform 40 year old ED Bays with “Minimal” Effort and Down Time

CHALLENGE:

- Hidden and Mismanaged Supplies
- Wasted Inventory
- Cluttered Floor Space
- Inflexible Wall-mounted Equipment
- Difficult to clean

When every second counts, it's critical to have what you need, where you need it, when you need it.

BEFORE



SOLUTION:

The “Evolution” equipment rail system

- Par level inventory that's quickly visible and easy to access
- Clear floor space that's easy to clean
- Equipment and supplies that can be moved at will for optimum workflow, without facility damage

AFTER




CHALLENGE:

Digging for supplies with gloved hands breeds issues of:

- Cross-contamination
- Hospital Acquired Infections
- Wasted time

Inventory systems are plagued by “First-In, Last-Out” syndrome. Unused materials become expired or cannot be restocked into Clean Supply creating wasted inventory costs.



SOLUTION:

Implementing the “Evolution” Medical Equipment Management Rail and Accessories System improves cost and work-flow efficiencies and infection control. It reduces the risks of cross-contamination, the acquisition of HAI’s, the wasted expense of expired and unused materials, and wasted time.

LEAN supply inventory management at the bedside provides clinicians with:

WHAT THEY NEED
WHEN THEY NEED IT
WHERE THEY NEED IT
IN QUANTITIES THEY NEED

Permanently reduce waste from over purchasing and from items that cannot be restocked into Clean Supply.

- Clips, Hooks, Bins, and Baskets of varying sizes and multi-use purposes
- Clear bins for quick visual identification of stock for supplies
- Embedded anti-microbial agents in all plastic adapters and powder coat painted finish
- Snap-on, snap-off features for easy relocation and reorganization of supplies
- Flexible, Adaptable, Expandable for additional organizational accessories as needed

CHALLENGE:

How can we ensure that we design and build the best possible solution for our clinical team partners?

Is this the best we can do?



Or this?




Of course not!
Yet it's the standard of practice.



SOLUTION:

Place all hand washing compliance and infection control accessories on the rail.

Studied and published by the Center for Health Design in 2012: the evidence-based recommendation is to create a standard directional flow with easy access to sink and hand washing accessories at the room entry by "...placing the hand washing accessories on an equipment rail allows adjustment in the system without creating damage to the wall and potential home for pathogens."



True Terminal Cleaning and Infection Control

CHALLENGE:

The dangers of contagious infectious disease are very costly problems that can have deadly results. With or without a rail system, if the walls are not stripped bare of all their elements, infectious diseases have places to hide. This method is unacceptable.




True Terminal Cleaning is only possible with the ability to remove ALL equipment and materials from the walls without leaving holes in the walls, without cleaning around wall mounted items, without leaving places for infection to hide.



If you can not remove it, you can not clean it!

SOLUTION:

Experts agree that the best Terminal Cleaning Solution requires bare walls, and free access to all surfaces.





True Terminal Cleaning looks like this.



Complete Design Services, including AutoCAD, REVIT / BIM and 3D Renderings!

**Call or email us today,
so we can take care of your new project!**



www.PaladinHC.com

Phone: 407 – 614 8781 - Email: info@PaladinHC.com

IF IT'S ON A WALL – IT'S ON A RAIL!



IF YOU NEED IT – WE CAN DO IT!

**Emergency
Department**

ICU

Patient Room

NICU

**Physician /
Clinic Exam Room**

Pre and Post Op



Vizient Supplier Number: CE7040
Premier Supplier Number: PP-FA-945

Scrub Sinks

**General Infection
Control**

Operating Rooms

Ambulatory Care

**Mass Casualty /
Surge Population
Management**

Forward Triage Cart