



**Mounting Plate**  
**with 1 hook hanger**  
**powder coated metal wire hook**

**Item Numbers:**  
**SM301-1-PM (1" hook)**  
**SM301-2-PM (2" hook)**  
**SM301-4-PM (4" hook)**

**Requires user defined  
JSN number**

**Dimensions:**  
**1.5"W x 2.5"H**  
**x length of hook"D, Weight: 1lb**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.