



Mounting Plate
with gloves box holder,
heavy duty ABS plastic with
internal dividers (glove boxes
are not included)

Available in single, double,
triple and quadruple. Please
specify prior to order
placement.

Item Number: IC227-PM
JSN Number: A5107A

Dimensions:
10”H x 18.75”W x 5”D (triple)
Weight: 3lbs

Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.