



**Mounting Plate**  
**with heavy duty acrylic gloves**  
**box holder**  
**(glove boxes are not included)**



**Available in single, double, triple and quadruple. Please specify prior to order placement.**

**Item Number: IC287-PM**  
**JSN Number: A5107A**

**Dimensions:**  
**10"W x 18.75"H x 5"D (triple)**  
**Weight: 3lbs**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.