



Rail Mounting Kit with IV swing-arm holder and 4-hook IV bag hanger, infusion pump mounting swing arm for attachment of infusion pump and other pole-mount devices. Includes adjustable telescoping IV pole with 4 bag hanger. Does not include GCX channel.

Item Numbers: IM301-K

Requires user defined JSN number

**Dimensions:
37"D x 5"H x 2"W
Weight: 12lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.