



**Mounting Plate**  
**for\_FERNO oxygen station**  
**#T70710-N with 2 oxygen tanks or**  
**without tanks (please specify at time**  
**of order), backwall drilled for Paladin**  
**adapter mounting**

**Item Numbers:**  
**FN508-P (Plate only)**  
**FN508-PM (Plate & Oxygen Station)**

**Requires user defined**  
**JSN number**

**Dimensions:**  
**18"W x 34.5"H x 10"D (Oxygen Station)**  
**Weight: 3lbs (Plate only)**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.