



**Mounting Plate**  
**with medical-gas block**  
**with 2 x 1/8" NPT male front outlets. Air and oxygen blocks come with an additional 90 degrees DISS center outlet check valve. Vacuum block does not include center check valve.**

**Item Numbers:**  
**GBA105 (Air)**  
**GBO105 (Oxygen)**  
**GBV105 (Vacuum)**  
**Requires user defined JSN number**

**Dimensions: 4"W x 3"H x 3.5"D**  
**Weight: 3lbs**

**Flow meters and vacuum regulators are not included**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.