



Medical Gas Block with 1 x 1/8" NPT Male fitting front outlet, 1 x Female Quick Connect front Outlet. Available with Ohio, Chemetron, Oxyquip and Schraeder fittings. Air and Oxygen blocks come with an additional 90 degrees DISS center outlet check valve. Vacuum block does not include center check valve.

Item Numbers:

GBA104 (Air)

GBO104 (Oxygen)

GBV104 (Vacuum)

JSN Number: not established

are not included

Dimensions: 4"W x 3"H x 3.5"D

Weight: 3lbs

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.