



**Rail Mounting Kit with**  
**VU Ease Patient Information Board**  
**The Information Boards are built**  
**with high-quality tempered glass and**  
**include a limited lifetime warranty.**

**Item Numbers:**

**AC149-14.75x12 (small)**

**AC149-16x26 (medium)**

**AC149-17.75x15 (large)**

**JSN Number: not established**

**Dimensions:**

**14.75"W x 12"L x 2"D (small)**

**16"W x 26"L x 2"D (medium)**

**17.75"W x 15"L x 2"D (large)**

**Weight: 6lbs**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.