



**Rail Mounting Kit**  
**with 2-tiers Medical Chart Holder**  
**Acrylic, Opaque front face,**  
**HIPAA compatible**

**Item Numbers:**  
**AC173-2-PM-HIPAA (2-tiers)**  
**JSN Number: not established**

**Dimensions:**  
**9.5”W x 14”H x 3.5”D**  
**Weight: 5lbs**

**Photo shows 3-tiers chart holder for reference**

**Item Description:**

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.