



**PART NUMBER:** MC101 (JSN # not established)

**DESCRIPTION:** F.T.S.  
Flexible Treatment Station

- NOTES:**
- Manage all clinical and support equipment in forward triage or hallway overflow conditions
  - Use as a flu screening cart in the ED waiting areas to alleviate throughput
  - Able to supply locking drawers and a variety of work surfaces and file organizers for admission forms, etc.
  - Stable weighted base that is powder-coated in an anti-microbial finish
  - Heavy-duty rubberized non-skid casters, 2 of which are lockable
  - (20) Linear Feet of usable equipment rail space for custom location of all accessories based on immediate need
  - Adjustable transport docking device will adapt to uneven surfaces or emergency Trendelenburg positioning
  - Includes (2) infusion poles and (5) EVOLUTION equipment rails
  - 26"W x 26"D x 60"H (without IV pole) – 77"H with max IV pole

