



**PART NUMBER:** MC102 (JSN # not established)

**DESCRIPTION:** Mass Casualty Rapid Deployment Equipment Management System

**NOTES:** Portable and easily stored Equipment Management System for rapid deployment in any situation. Application on rough terrain is accommodated through use of heavy duty Tri-Pod Stand which will “self level”.

- Includes (2) tri-pods with (2) carrying cases (2) 4ft rails with rail joiner in canvas pouch
- Sturdy and stable
- Pre-Hospital or In-Hospital
- Perfect for MCI Tent Applications
- Can be stored in a closet, MCI Trailer or heavy rescue type vehicle
- Delivers Med Gas, mounts PPE and monitoring devices on rail and off the floor
- Complete Integration with FERNO Intraxx components

