



PART NUMBER: SM501-4-K (Requires user defined JSN number)

DESCRIPTION: Small Cart

NOTES:

- Manage all clinical and support equipment in forward triage or hallway overflow conditions
- Use as a flu screening cart in the ED waiting areas to alleviate throughput
- Able to supply locking drawers and a variety of work surfaces and file organizers for admission forms, etc.
- Stable weighted base that is powder-coated in an anti-microbial finish
- Heavy-duty rubberized non-skid casters, 2 of which are lockable
- (20) Linear Feet of usable equipment rail space for custom location of all accessories based on immediate need
- Adjustable transport docking device will adapt to uneven surfaces or emergency Trendelenburg positioning
- Includes (2) infusion poles and (5) EVOLUTION equipment rails
- 26"W x 26"D x 60"H (without IV pole) – 77"H with max IV pole

