



## Mounting Plate with Cord Caddy

**Cord / hose storage and retrieval system**



**Item Number:**

**AC126-PM-2 (double)**

(8.5”H x 5.5”W x 2.5”D)

**AC126-PM-3 (triple)**

(8.5”H x 9”W x 2.5”D)

**AC126-PM-4 (quadruple)**

(8.5”H x 12”W x 2.5”D)

**AC126-PM-5 (quintuple)**

(8.5”H x 15”W x 2.5”D)

**AC126-PM-6 (sextuple)**

(8.5”H x 18”W x 2.5”D)

**Weight 3lbs**

**Requires user defined**

**JSN number**

### Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.