



**GCX101-25-2 Channel to span
(2) rails at 12-22.5" between rails**
**GXC101-31-2 channel to span
(2) rails at 15-28" between rails**
Additional lengths available,
To mount a variety of equipment,
arms, and management tools.
Includes (2) pre-installed
adapters.

Item Number:

GCX101-25-2 (25"H x 4"W x 1"D)

Weight 10lbs

GCX101-31-2 (31"H x 4"W x 1"D)

Weight 12lbs

JSN Number: not established

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.