



**Rail Mounting Kit with forward or backwards facing cable / gas hose management adapter**  
**Accommodates cables and gas hoses for secure and orderly management. Powder coated antimicrobial finish (AMF)**

**Item Number:**

**SM311 – forward facing**

**SM311 – backward facing**

**JSN Number: not established**

**Dimensions:**

**1.5”W x 1.25”D x 3.5”H**

**Weight: 1lb**

**Item Description:**

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.