



**Universal Acrylic Face Mask
with Face Shield Dispenser
Holder, including rail mounting
kit, heavy duty acrylic**

**For a complete overview which
brand items will fit, please see
next page or visit
www.PaladinHC.com**



**Item Number: IC-PFM200-PM
JSN Number: not established**

**Dimensions:
9.625"W x 8.5"H x 4.625"D
Weight: 3lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.



PART NUMBER: IC-PFM200-PM

DESCRIPTION: Universal Acrylic Face Mask with Face
Shield Dispenser Holder

NOTES: 9.625"W x 4.625"D x 8.5"H (dispenser only)
Heavy Duty Acrylic, Dispenser only,
incl. Rail mount Kit.

**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1120/37062, 1832, 1838R, 1870, 1870-N95, 8214, 8514, 9210/37021

Cardinal Health: AT7505-P, AT7511-WE, AT51035, AT52005, AT54635,
AT71035, AT71039, AT72835, AT74535, AT74631, AT74635, AT771145,
N95-ML, 22EL77

Halyard: 28804, 28809, 28821, 37505, 37510, 41803, 47095, 47137,
47147, 47148, 47149, 47298, 47500, 48225, 48227, 48228, 48237, 48247,
48248, 48257, 48298, 49210, 49212, 49230, 59926, 59930, 59937, 62114,
62116, 62368, 69801

McKesson: 91-1600, 91-2600

MediChoice, Primagard/PriMed: 41573, 41592, 42073, 42333, 42560,
42562, 43092, 43573, 44073, 45033, 45073, 46003, 46403, 77701, 77702,
77703, PG43573

Medline: NON27373A, NON27410, NON27410EL, NON27411,
NON27411AG, NON27600, NON27600AG, NON27710, NON27710EL,
NON27810, NON27810EL, NON27361A, NON27382, NON27408,
NON27712

Precept /Prestige Ameritech/ProGear: 65-3344, 15310, 15330, 15710

Sperian: 091007, 14110452