



Universal Face Mask with Ear loop Dispenser Holder, including rail mounting kit, heavy duty acrylic

For a complete overview which brand items will fit, please see next page or visit www.PaladinHC.com

**Item Number: IC-PFP022-PM
JSN Number: not established**

**Dimensions:
8.125"W x 7.5"H x 5.625"D
Weight: 3lbs**

Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.



PART NUMBER: IC-PFP022-PM

DESCRIPTION: Universal Face Mask with
Earloops Dispenser Holder

NOTES: 8.125"W x 5.625"D x 7.5"H (dispenser only)
Heavy Duty Acrylic, Dispenser only,
incl. Rail Mount Kit.

**FITS THE
FOLLOWING
DISPENSERS:**

3M: 1832, 1838, 1900, 8214, 8514

Cardinal Health: AT7505-P, AT51035, AT71035, AT72835, AT74535,
AT752005, AT771145, N95-ML

Halyard: 28800, 28804, 28809, 28820, 28821, 37505, 37510, 41803,
47137, 47147, 47148, 47298, 47500, 47625, 47650, 48207, 48208,
48210, 48220, 48225, 48227, 48228, 48237, 48247, 48248, 48257,
48297, 48298, 48390, 49210, 49212, 49215, 49216, 49230, 49235,
59926, 59930, 59937, 62114, 62116, 62368

McKesson: 91-1600, 91-2600

MediChoice, Primagard/PriMed: PG42073, PG42092, PG42121,
PG42332, PG42333, PG43573

Medline: NON27361A, NON27371A, NON27378A, NON27379A,
NON27382, NON27382FR, NON27408, NON27420, NON27420EL,
NON27712

Precept: 15212, 15216, 15220, 15310, 15330, 15710, 65-3344

Prestige Ameritech/ProGear: FM81110