




**Universal Face Mask with
Face Shield Dispenser Holder,
including rail mounting kit,
heavy duty acrylic**

**For a complete overview which
brand items will fit, please see
next page or visit
www.PaladinHC.com**



**Item Number: IC-PFP120-PM
JSN Number: not established**

**Dimensions:
9.125"W x 7.5"H x 4.25"D
Weight: 3lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.



PART NUMBER: IC-PFP120-PM

DESCRIPTION: **Universal Face Mask with Face Shield
Dispenser Holder**

NOTES: 9.125"W x 4.25"D x 7.5"H (dispenser only)
Heavy Duty Acrylic, Dispenser only,
incl. Rail Mount Kit.

**FITS THE
FOLLOWING
DISPENSERS:**

Halyard:
28804, 28809, 28821, 41803, 47137, 47147, 47148, 47149,
47298, 47500, 48228, 48237, 48247, 48248, 48257, 48298,
59930, 59937, 62114, 62116, 69801

McKesson:
91-1600, 91-2600

MediChoice, Primagard/PriMed:
PG41573, PG41592, PG45033, PG45073

Medline:
NON27600, NON27600AG, NON27710, NON27710EL, NON27810,
NON27810EL

Precept:
65-3344, 15310, 15330, 15710

Prestige Ameritech/ProGear:
FM82120, FM84210

Sperian:
14110452