



**Universal Acrylic Face Mask  
with Ear loop Dispenser Holder,  
including rail mounting kit,  
heavy duty acrylic**

**For a complete overview which  
brand items will fit, please see  
next page or visit  
[www.PaladinHC.com](http://www.PaladinHC.com)**



**Item Number: IC219-PM  
JSN Number: not established**

**Dimensions:  
8.125"W x 4.75"H x 5.625"D  
Weight: 3lbs**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.



---

**PART NUMBER:** IC219-PM

**DESCRIPTION:** **Universal Acrylic Face Mask  
with Earloop Dispenser Holder**

**NOTES:** 8.125"W x 5.625"D x 4.75"H (dispenser only)  
Heavy Duty Acrylic, Dispenser only,  
incl. Rail Mount Kit.

---

**FITS THE  
FOLLOWING  
DISPENSERS:**

**3M:** 1800, 1816, 1820, 1826, 1832, 1838, 1860S

**Cardinal Health:** AT70021, AT71021, AT771141, AT771211

**Halyard:** 25867, 25868, 25869, 28797, 28802, 28806, 28808, 28810,  
32856, 37525, 37525, 39117, 41802, 41805, 47080, 47085, 47090,  
47107, 47108, 47117, 47119, 47127, 47295, 47297, 48100, 48105,  
48201, 48215, 48296, 48310, 49214, 49310, 49700, 49701, 59927,  
59928, 59929, 59931, 59936, 59938, 59939, 59940, 62113, 62115,  
62356, 62357, 62362, 62363, 62364, 62367, 62369, V1025,

**McKesson:** 91-1000, 91-1200, 91-1400, 91-2002

**MediChoice, Primagard/PriMed:** 1773ELM001, 41001, 41073, 41092,  
41200, 41201, 41241, 41331, 42000, 42001, 42092, 42121, 42241,  
42320, 42331, 42332, 42421, 44021, 44051, 44241, 44331, 46001

**Medline:** NON27110, NON27120, NON27120-2, NON27122,  
NON27355, NON27365, NON27374, NON27375, NON27376,  
NON27377, NON27378, NON27385, NON27386, NON27400,  
NON27402, NON27408EL, NON27412, NON27412EL, NON27712EL,  
NON27812EL, NONSH400

**Precept Prestige Ameritech/ProGear:** 15100, 15101, 15130, 15150,  
15155, 15160, 15200, 15211, 15212, 15230, 15240, 15260, 15300,  
15320, 15420, 15500, 65-3320, 65-3322, FM81000, FM85011

**Tidi:** 9010, 9020, 9030, 9040