



*Customer to define mounting preference
(horizontal or vertical) at time of order.*

**Rail Mounting Kit with PPE
Dispensing Cabinet, horizontal
top feed application, ABS
plastic, holds 2 boxes of gloves,
1 box of face masks and 1
bottle of hand sanitizer**

**Item Number: IC271-PM
JSN Number: A5107**

**Dimensions:
18.5"W x 11.625"H x 4.5"D
Weight: 6lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.