



*Customer to define mounting preference (horizontal or vertical) at time of order.*

**Rail Mounting Kit with PPE Dispensing Cabinet, horizontal top feed application, ABS and PETG plastic, holds 3 boxes of gloves, 2 boxes of face masks, multi fold towels and patient bibs. Taps on both sides to hang patient bib clips.**

**Item Number: IC274-PM  
JSN Number: A5107**

**Dimensions:  
17.625"W x 23.5"H x 5"D  
Weight: 7lbs**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.