



**Rail Mounting Kit,**  
**for owner provided**  
**GA Pacific SofPull paper towel**  
**dispensers (models 52109,**  
**53009, 58201, 58217 and**  
**58008)**

**Specific model number must be**  
**verified prior to order**  
**placement.**

**Item Number: IC408-P**  
**JSN Number: A5080-RM**

**Dimensions:**  
**4.5"W x 0.125"D x 6.5"H**  
**Weight: 2lbs**

**\*\* Does not include dispenser A5080**

**Item Description:**

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.