



**Rail Mounting Kit for**  
**Owner provided Welch Allyn**  
**71420 tray for**  
**Welch Allyn 71140 Desk**  
**charger**

**Item Number: DEM103-P**  
**JSN Number: not established**

**Dimensions:**  
**1.5”H x 7.5”W x 4.75”D**  
**Weight: 4lbs**

Adapter plate includes pre-installed adapters  
and wall stand-offs

**\*\* Welch Allyn equipment SOLD SEPARATELY**

**Item Description:**

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.