



**Adapter**  
**for owner provided Exergen**  
**Temporal Scanner Thermometer**  
**TAT5000 or wall mount Exergen**  
**134201**  
**(Exergen products are not**  
**included)**

**Item Number: DEM146-A**

**Requires user defined**  
**JSN number**

**Dimensions:**  
**2.6”H x 3”W x 1”D**  
**Weight: 0.5lb**

**Item Description:**

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.