



MADE IN THE USA



** Welch Allyn
equipment SOLD SEPARATELY

Item Description:

- Mounts to “EVOLUTION,” Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer’s devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.

Mounting Plate for owner provided Welch Allyn S- shaped spot monitors with horizontal monitor mount (Welch Allyn products are not included)

Welch Allyn Spot Vital Signs Device with Nellcor SpO2, SureTemp Thermometry, Non-invasive Blood Pressure, Pulse Rate, MAP; 120 V, 60 Hz AC, Lithium Ion (#105631) Battery; IEC Plug Type-B. Various part numbers available for this monitor type. Part Number: 42NTB-E1 - GTIN: 00732094005714

Item Number: DEM164-P
JSN Number: M4116-RM

Dimensions:
12”H x 3”W x 1”D
Weight: 1lb