



Mounting Plate with fixed shelf

with lip edges, metal with antimicrobial (AMF) powder coat painted finish, includes 0.75" lip on 3 sides (Monitor not included)

Item Number:

S201-PM (small)

S202-PM (medium)

S203-PM (large)

Requires user defined JSN number

Dimensions:

11"W x 8.5"D (small shelf overall outside dimensions)

10"W x 8"D (small shelf usable space)

11"W x 13.5"D (medium shelf overall outside dimensions)

10"W x 13"D (medium shelf usable space)

14"W x 13.5"D (large shelf total dimensions)

13" W x 13"D (large shelf usable space large shelf)

Weight: 3lbs (small) / 4lbs (medium) / 5lbs (large)

Weight capacity: 5lbs



Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.