



Rail Mounting Kit with fixed shelf with lip edges, factory attached, Metal with antimicrobial (AMF) powder coat painted finish, includes 0.75" lip on 3 sides

Item Number:

S201 (small)

S203 (large)

JSN Number: not established

Dimensions:

10"W x 8"D (small shelf)

14.5" W x 13.5"D (large shelf)

Weight: 3lbs (small) / 4lbs (large)

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.