



Mounting Plates
with heavy duty flip shelf,
scratch resistant food grade high density polyethylene. Intended spill lip. Up to 100lbs load capacity. (Item on shelf for illustration purposes only. Not included.)



Item Number:
S242-PM (small shelf)
S241-PM (medium shelf)
S244-PM (large shelf)

Requires user defined JSN number

Dimensions:
15"W x 10"D – 4lbs (small shelf)
20"W x 13"D – 5lbs (medium shelf)
24"W x 18"D – 6lbs (large shelf)

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.