



Mounting Plate with holder
for owner provided round
bottle type dispensers / hand
pumps 2.5" and 3.25" diameter

e.g. 3M Avagard Nail Cleaner
#9204, Aero Wash Eye Rinse
#AW160ZEPUS, Cederroth
#725200A, Pure Nest Dispenser

Item Number: IC141-PM

**Requires user defined
JSN number**

**Dimensions:
6.4"H x 3.3"W x 2.9" / 3.25"D
Weight: 2lbs**

Item Description:

- Mounts to "EVOLUTION," Fairfield-style Medical Equipment Management Rail using a snap action adapter with positive snap attachment lock feature.
- Integrated Allen screws provide the ability to lock the adapter in place on rail.
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.